

Product datasheet

Anti-beta Defensin 1 antibody ab115813

1 References 1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-beta Defensin 1抗体 |
| 描述 | 兔多克隆抗体to beta Defensin 1 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB, IHC-P |
| 种属反应性 | 与反应: Human |
| 免疫原 | Synthetic peptide corresponding to a sequence at the C terminal of Human beta Defensin 1. |
| 阳性对照 | Human lung cancer tissue. |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservatives: 0.025% Thimerosal (merthiolate), 0.025% Sodium azide Constituents: 2.5% BSA, 0.45% Sodium chloride, 0.1% Dibasic monohydrogen sodium phosphate |
| 纯度 | Protein A purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab115813** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab 评论 | 说明 |
|----|-------|----|
|----|-------|----|

| | |
|----|--|
| WB | Use a concentration of 1 µg/ml. Predicted molecular weight: 7 kDa. The detection limit for beta Defensin 1 is approximately 1ng/lane under non-reducing and reducing conditions. |
|----|--|


| 应用 | Ab 评论 | 说明 |
|----|-------|----|
|----|-------|----|

IHC-P Use a concentration of 1 µg/ml. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

靶标

| | |
|-------|--------------------------------------|
| 功能 | Has bactericidal activity. |
| 组织特异性 | Plasma. |
| 序列相似性 | Belongs to the beta-defensin family. |
| 细胞定位 | Secreted. |

图片



ab115813 at 1µg/ml staining beta Defensin 1 in by Paraffin-embedded Human Lung cancer tissue by Immunohistochemistry,

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-beta Defensin 1 antibody (ab115813)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors