

Product datasheet

Anti-CIP75 antibody ab106443

★★★★☆ 2 Abreviews 2 References 2 图像

概述

产品名称	Anti-CIP75抗体
描述	兔多克隆抗体to CIP75
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide corresponding to Human CIP75 (N terminal). Database link: Q9NRR5
阳性对照	3T3 cell lysate, Mouse kidney tissue.
常规说明	Previously labelled as UBQLN4.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
存储溶液	Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106443** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml. Perform heat mediated antigen retrieval via the microwave method before commencing with IHC staining protocol.


应用	Ab评论	说明
----	------	----

WB ★★★★★ Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 64 kDa.

靶标

组织特异性	Highly expressed in pancreas, kidney, skeletal muscle, heart and throughout the brain, and at lower levels in placenta, lung and liver.
序列相似性	Contains 1 UBA domain. Contains 1 ubiquitin-like domain.
翻译后修饰	Phosphorylated upon DNA damage, probably by ATM or ATR.
细胞定位	Nucleus. Cytoplasm. Endoplasmic reticulum. May associate with the endoplasmic reticulum.

图片



Western blot - Anti-CIP75 antibody (ab106443)


Lane 1 : Anti-CIP75 antibody (ab106443) at 1 µg/ml

Lane 2 : Anti-CIP75 antibody (ab106443) at 2 µg/ml

All lanes : 3T3 cell lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 64 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CIP75 antibody (ab106443)

ab106443, at 5 µg/ml, staining CIP75 in Mouse kidney tissue by immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors