

Product datasheet

FITC Anti-6X His tag® antibody ab1270

1 References

概述

产品名称	FITC荧光Anti-6X His tag®抗体
描述	FITC荧光山羊多克隆抗体to 6X His tag®
宿主	Goat
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: ICC
种属反应性	与反应: Species independent
免疫原	Synthetic peptide: HHHHHH (6-His) conjugated to KLH.

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

HIS-TAG® is a trademark of EMD Biosciences, Inc.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium azide Constituent: 0.2% BSA
纯度	Immunogen affinity purified
纯化说明	Antibodies were affinity purified using the peptide immobilized on solid support. Antibody concentration was determined by extinction coefficient: absorbance at 280 nm of 1.4 equals 1.0 mg of IgG.
克隆	多克隆

应用

The Abpromise guarantee [Abpromise™](#) 承诺保证使用ab1270于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC		1/200 - 1/2000.

靶标

相关性 The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. His-tags can be fused to the amino- or carboxy- termini of proteins in transfected or transformed cells.

细胞定位 Depends upon the localization of the parent protein tagged with hexahistidine.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors