

Product datasheet

FITC Anti-CD5 antibody [CRIS1] ab1156

1 References

概述

产品名称	FITC荧光Anti-CD5抗体[CRIS1]
描述	FITC荧光小鼠单克隆抗体[CRIS1] to CD5
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to CD5.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.4 Preservative: 0.1% Sodium azide Constituents: 0.5% BSA, PBS
克隆	单克隆
克隆编号	CRIS1
骨髓瘤	unknown
同种型	IgG2a
轻链类型	unknown

应用

应用说明	1. Characterization of leukaemias and lymphomas in human lysed whole peripheral blood or
-------------	--

mononuclear cells separated by density gradient.

2. In combination with CD19 (ab), this antibody may be used for the diagnosis of B-chronic lymphoblastic leukaemia.

CD5 (FITC) immunofluorescence analysis can be performed on a flow cytometer equipped with an excitation source of 488nm and fitted with logarithmic amplifiers.

10µl of CD5 (FITC) is sufficient for labelling of 1×10^6 cells.

Molecular F/P ratio: between 3.0 and 6.0

靶标

功能	May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.
序列相似性	Contains 3 SRCR domains.
细胞定位	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors