

Anti-Hepatitis E Virus antibody ab1036

概述

产品名称	Anti-Hepatitis E Virus抗体
描述	兔多克隆抗体to Hepatitis E Virus
宿主	Rabbit
特异性	This antibody recognizes the capsid protein of the HEV variant HEV US-u3-2e identified in the United States.
经测试应用	适用于: WB
种属反应性	与反应: Other species
免疫原	Peptide corresponding to amino acids 627-637 (referred to as US-u3-2e)of the 71 kD Hepatitis E Virus (HEV) capsid protein (US-1 strain of HEV). This peptide sequence is highly conserved in human and swine and in the Burmese, Mexican, and Japanese strains of HEV.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.4 Preservative: 0.05% Sodium azide
纯度	Immunogen affinity purified
克隆	多克隆
同种型	unknown
轻链类型	unknown

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab1036于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000.

靶标

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors