

Product datasheet

Mouse IgG1, kappa monoclonal [15-6E10A7] - Isotype Control
ab170190

17 References 3 图像

概述

产品名称	小鼠IgG1, kappa单克隆抗体[15-6E10A7]-同型对照
特异性	This Mouse IgG1 Isotype Control Antibody was raised in mouse against a yeast-specific protein making it unsuitable for yeast experiments.
经测试应用	适用于: Flow Cyt, ICC, IHC-P
免疫原	Recombinant full length protein. ab170190 binds specifically to DNA-binding domain of the yeast GAL4 protein transcription activator.
常规说明	<p>This antibody clone is manufactured by Abcam.</p> <p>Isotype controls are used to confirm that the primary antibody binding is specific and not a result of non-specific Fc receptor binding or other protein interactions. The isotype control antibody should match the primary antibody's host species, isotype, and possible conjugation. The control performed appropriately in all materials and platforms that were tested.</p> <p>Product was previously marketed under the MitoSciences sub-brand.</p> <p>If you require this antibody in a particular buffer formulation or a particular conjugate for your experiments, please contact orders@abcam.com or you can find further information here.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.5</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: HEPES buffered saline</p>

纯度	Ammonium Sulphate Precipitation
克隆	单克隆
克隆编号	15-6E10A7
同种型	IgG1
轻链类型	kappa


应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab170190于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Flow Cyt		Use a concentration of 1 µg/ml.
ICC		Use a concentration of 1 µg/ml.
IHC-P		Use at an assay dependent concentration. PubMed: 26627833

图片



Flow cytometry experiments with 4% PFA fixed HeLa (Human adenocarcinoma), 653s (mouse myeloma), H4IIE (rat hepatoma) and H9C2 (rat myoblast) were performed with the IgG1 Isotype Control Antibody (red) and no primary antibody negative control (black).


An Alexa Fluor® 488 conjugate with isotype specificity to the mouse antibody was used as a secondary antibody. The isotype control at 1 ug/mL shows no higher signal than the no primary negative control.



Immunocytochemistry - Mouse IgG1, kappa monoclonal [15-6E10A7] - Isotype Control (ab170190)

Immunocytochemistry (ICC) experiments with NIH-3T3 (mouse sarcoma), COS7 (monkey kidney fibroblast), H9C2 (rat myoblast) and HeLa (human adenocarcinoma) were performed with the IgG1 Isotype Control Antibody (top), no primary antibody negative control (middle), and [ab14705](#) as a positive control (bottom).

An Alexa Fluor® 488 conjugate with isotype specificity to the mouse antibody was used as a secondary antibody. The isotype control at 1 µg/mL shows no higher signal than the no primary negative control.



ELISA - Mouse IgG1, kappa monoclonal [15-6E10A7] - Isotype Control (ab170190)

Fluorescence ELISA.

An isotyping ELISA was performed by coating a 96-well plate with 1.25 µg/mL of the IgG1 Isotype Control Antibody and detecting with Alexa Fluor conjugates specific to mouse IgG1, IgG2a, IgG2b, IgG3, IgM and heavy and light chains (H&L) of IgG. This experiment verifies that the primary antibody's isotype is correct and that it is successfully bound by the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors