

Product datasheet

Anti-p66 beta antibody ab13713

1 图像

概述

产品名称	Anti-p66 beta抗体
描述	兔多克隆抗体to p66 beta
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide, corresponding to amino acids 445-460 of Mouse p66 beta .
阳性对照	Raw or HeLa whole cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS, 0.1mg/ml BSA
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab13713** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 66 kDa.

靶标

功能 Has transcriptional repressor activity. Targets MBD3 to discrete loci in the nucleus.

组织特异性	Ubiquitous.
序列相似性	Contains 1 GATA-type zinc finger.
结构域	CR1 is required for interaction with MBD3 and transcription repression. CR2 is required for localization to discrete loci in the nucleus.
细胞定位	Nucleus. Nuclear, in discrete foci.

图片



Western blot analysis of p66b in 15 ugs of Raw cell lysate using ab13713 at 1/500 dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors