



Product datasheet

Normal Mouse Serum (sterile) ab7486

13 References

概述

产品名称	Normal小鼠Serum (sterile)
宿主	Mouse
经测试应用	适用于: Blocking, IHC-P, IHC-Fr
常规说明	Mouse serum for blocking in antibody-based applications or use as an antibody diluent. Origin: North America. Total Protein Concentration : normal Strain: Mouse - Mixed Sex: Mixed sex Raised in: Mouse Purity: Whole antiserum Sterilization: This product was aseptically filtered through a 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.

性能

形式	Liquid
存放说明	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: Whole serum

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7486于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Blocking		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors