



Product datasheet

Normal Sheep Serum (sterile) ab7489

2 References

概述

产品名称	Normal羊Serum (sterile)
宿主	Sheep
经测试应用	适用于: IHC-Fr, Blocking, IHC-P
常规说明	<p>Sheep serum for blocking in antibody-based applications or use as an antibody diluent. From mixed sex mixed strain sheep. Origin: North America.</p> <p>Total Protein Concentration : 90 mg/ml (by Refractometry). Sex: Mixed breed and sex.</p> <p>Raised in: Sheep</p> <p>Purity: Whole antiserum</p> <p>Sterilization: This product was aseptically filtered through a 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.</p>

性能

形式	Liquid
存放说明	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab7489于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr		Use at an assay dependent concentration.
Blocking		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors