

10x EDTA Buffer pH 8.0 ab64216

2 References 1 图像

概述

|       |  |
|-------|--|
| 产品名称  | 10x EDTA缓冲液pH 8.0  |
| 经测试应用 | 适用于: IHC-P   |
| 常规说明  | EDTA buffer pH 8.0 for heat-induced antigen retrieval (HIER) during IHC (ab64216).<br><br><b>Recommended Procedure: Microwave</b><br><br>Deparaffinize and rehydrate tissue sections.<br>Dilute the product 1:10 with deionized water.<br>Immerse sections in coplin jar filled with target retrieval solution microwave on HIGH power until boiling.<br>Keep the slides warm by heating for 10min at LOW power.<br>Let coplin jar sit in microwave for at least 20min.<br>Remove and rinse with buffer.<br>Proceed with staining. |

性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Store at +4°C.   |
| 存储溶液 | pH: 7.7<br>Constituent: 0.29% EDTA<br><br>Buffer 10x concentrated.   |
| 相关性  | EDTA buffer is used on formalin-fixed, paraffin-embedded tissue sections mounted on glass slides for target retrieval prior to immunohistochemistry (IHC) procedures. Target retrieval prior to IHC procedures obtains positive results, or results in an increase in staining intensity with many primary antibodies. |


应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab64216于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用    | Ab评论 | 说明                                  |
|-------|------|-------------------------------------|
| IHC-P |      | Use at an assay dependent dilution. |

## 图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - 10x EDTA Buffer pH 8.0  
(ab64216)

ab64216 Human Melanoma stained with an anti-S-100 antibody diluted with Antibody Diluent ([ab64211](#))

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors