

Bag1 overexpression 293T lysate (whole cell) ab94217

2 图像

概述

产品名称	Bag1 overexpression 293T裂解物(whole cell)
常规说明	ab94217 is a 293T cell transfected lysate in which Human Bag1 has been transiently over-expressed using a pCMV-Bag1 plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Inhibits the pro-apoptotic function of PPP1R15A, and has anti-apoptotic activity. Markedly increases the anti-cell death function of BCL2 induced by various stimuli. Tissue specificity: Isoform 4 is the most abundantly expressed isoform. It is ubiquitously expressed throughout most tissues, except the liver, colon, breast and uterine myometrium. Isoform 1 is expressed in the ovary and testis. Isoform 4 is expressed in several types of tumor cell lines, and at consistently high levels in leukemia and lymphoma cell lines. Isoform 1 is expressed in the prostate, breast and leukemia cell lines. Isoform 3 is the least abundant isoform in tumor cell lines (at protein level). Similarity: Contains 1 BAG domain. Contains 1 ubiquitin-like domain. PTM: Ubiquitinated; mediated by SIAH1 or SIAH2 and leading to its subsequent proteasomal degradation.


应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab94217于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent dilution.

图片



SDS-PAGE - Bag1 overexpression 293T lysate
(whole cell) (ab94217)

ab94217 at 15 μ g/lane on an SDS-PAGE gel.



Western blot - Bag1 overexpression 293T lysate
(whole cell) (ab94217)

All lanes : Anti-Bag1 antibody (**ab54817**) at 1/500 dilution

Lane 1 : Bag1 293T Transfected Lysate - Positive Control
(ab94217)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 μ g per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugate at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors