

### Cathepsin G overexpression 293T lysate (whole cell) ab94086

#### 2 图像

#### 概述

产品名称	Cathepsin G overexpression 293T裂解物(whole cell)
常规说明	ab94086 is a 293T cell transfected lysate in which Human Cathepsin G has been transiently over-expressed using a pCMV-Cathepsin G plasmid. The lysate is provided in 1X Sample Buffer.
经测试应用	适用于: WB

#### 性能

Mycoplasma free	Yes
形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Serine protease with trypsin- and chymotrypsin-like specificity. Cleaves complement C3. Has antibacterial activity against the Gram-negative bacterium P.aeruginosa, antibacterial activity is inhibited by LPS from P.aeruginosa, Z-Gly-Leu-Phe-CH <sub>2</sub> Cl and phenylmethylsulfonyl fluoride. Similarity: Belongs to the peptidase S1 family. Contains 1 peptidase S1 domain.

#### 应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab94086于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use at an assay dependent dilution.

#### 图片



SDS-PAGE - Cathepsin G overexpression 293T lysate (whole cell) (ab94086)

ab94086 at 15µg/lane on an SDS-PAGE gel.



Western blot - Cathepsin G overexpression 293T lysate (whole cell) (ab94086)

**All lanes** : Anti-Cathepsin G antibody (**ab89593**) at 1/500 dilution

**Lane 1** : Cathepsin G overexpression 293T lysate (whole cell) (ab94086)

**Lane 2** : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

**Secondary**

**All lanes** : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors