

Alkaline Phosphatase chromogen (BCIP/TNBT) ab7413

概述

产品名称	Alkaline Phosphatase chromogen (BCIP/TNBT)
经测试应用	适用于: WB, IHC-P, IHC-Fr
常规说明	<p>This product is a liquid substrate for alkaline phosphatase (AP), and is especially useful for the detection of AP-conjugated antibodies (i.e. secondary antibodies). Apply to tissue section for 10 minutes or until a suitable level of colouration has been achieved.</p> <p>BCIP: 5-bromo-4-chloro-3-indolyl phosphate 0.452 mM TNBT: tetranitroblue tetrazolium 0.432 mM This product used to be called "Alkaline phosphatase detection kit (BCIP/TNBT)".</p> <p>BCIP/TNBT is similar to BCIP/NBT but is more sensitive and produces a more intense purple color. A useful secondary marker in multiple staining procedures.</p>

性能

形式	Liquid
存放说明	Store at +4°C.
存储溶液	Constituents: 0.02% BCIP, 0.04% TNBT
Reagent说明	BCIP/TNBT is similar to BCIP/NBT but is more sensitive and produces a more intense purple color. A useful secondary marker in multiple staining procedures.

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7413于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors