

Anti-ABCA5 antibody ab99953

[2 References](#) [2 图像](#)

概述

产品名称	Anti-ABCA5抗体
描述	兔多克隆抗体to ABCA5
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant fragment corresponding to Human ABCA5 aa 1450 to the C-terminus. NM_018672 Database link: Q8WWZ7-1
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Lyophilized:Reconstitute with 200ul distilled sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 1% BSA
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab99953于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 187 kDa.
IHC-P		1/100 - 1/500.

靶标

功能	May play a role in the processing of autolysosomes.
组织特异性	Ubiquitously expressed. Highly expressed in testis, skeletal muscle, kidney, liver and placenta.
序列相似性	Belongs to the ABC transporter superfamily. ABCA family. Contains 2 ABC transporter domains.
发展阶段	Expressed in fetal liver, kidney and brain.
翻译后修饰	N-glycosylated.
细胞定位	Lysosome membrane. Late endosome membrane. Golgi apparatus membrane.


图片



Western blot - Anti-ABCA5 antibody (ab99953)

Anti-ABCA5 antibody (ab99953) at 1/500 dilution + Human fetal testis lysate

Predicted band size: 187 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ABCA5 antibody (ab99953)

ab99953 at 1/100 dilution, staining ABCA5 in human fetal kidney tissue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors