

Alkaline Phosphatase Anti-Biotin antibody ab106217

概述

产品名称	Alkaline Phosphatase Anti-Biotin抗体
描述	Alkaline Phosphatase兔多克隆抗体to生物素
宿主	Rabbit
偶联物	Alkaline Phosphatase
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase (calf intestine), anti-Rabbit Serum and Biotinylated IgG.
经测试应用	适用于: WB, ELISA, IHC-P
种属反应性	与反应: Species independent
免疫原	KLH conjugated Biotin
常规说明	<p>Dilute only prior to immediate use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Avoid freeze / thaw cycle.
存储溶液	<p>Preservative: 0.1% Sodium azide</p> <p>Constituents: 0.88% Sodium chloride, 1% BSA, 0.01% Magnesium chloride, 0.001% Zinc chloride, 50% Glycerol (glycerin, glycerine), 0.79% Tris HCl</p>
纯度	Immunogen affinity purified
纯化说明	ab106217 was purified from monospecific antiserum by immunoaffinity chromatography using Biotin coupled to sepharose beads.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab106217于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/2500.
ELISA		1/20000.
IHC-P		1/200 - 1/1000.

靶标

相关性

Biotin is a water soluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the initial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is required by all organisms but can only be synthesized by bacteria, yeasts, molds, algae, and some plant species. Biotin is required as prosthetic group of enzymes involved in incorporation of carbon dioxide into organic compounds. Biotin has a MW of 244 Da.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors