

Anti-Alpha Dystroglycan antibody [2238] ab106110

5 References [1 图像](#)

概述

产品名称	Anti-Alpha Dystroglycan抗体[2238]
描述	小鼠单克隆抗体[2238] to Alpha Dystroglycan
宿主	Mouse
特异性	ab106110 is specific for a glycoepitope present on brain Alpha Dystroglycan. The clone was generated using enriched bovine brain Alpha dystroglycan, and has been shown to bind Alpha Dystroglycan from enriched mouse, rat, rabbit and bovine brain tissues but did not bind Alpha Dystroglycan from striated muscle, peripheral nerve or kidney (PMID: 16709410).
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Rat, Human
免疫原	Full length native protein (purified) corresponding to Cow Alpha Dystroglycan. Purified bovine brain alpha-dystroglycan Database link: O18738
表位	Glycoepitope on brain bovine alpha-dystroglycan
阳性对照	WB: Human brain tissue lysate; SH-SY5Y cell lysate; Mouse brain tissue lysate; Rat brain tissue lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituents: 0.1% BSA, PBS
纯度	Protein G purified

纯化说明	0.2 μm filtered antibody solution
克隆	单克隆
克隆编号	2238
同种型	IgG2b

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab106110于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/100. Predicted molecular weight: 97 kDa.

靶标

功能	<p>The dystroglycan complex is involved in a number of processes including laminin and basement membrane assembly, sarcolemmal stability, cell survival, peripheral nerve myelination, nodal structure, cell migration, and epithelial polarization.</p> <p>Alpha-dystroglycan is an extracellular peripheral glycoprotein that acts as a receptor for both extracellular matrix proteins containing laminin-G domains, and for certain adenoviruses. Receptor for laminin-2 (LAMA2) and agrin in peripheral nerve Schwann cells. Also acts as a receptor for M.leprae in peripheral nerve Schwann cells but only in the presence of the G-domain of LAMA2, and for lymphocytic choriomeningitis virus, Old World Lassa fever virus, and clade C New World arenaviruses.</p> <p>Beta-dystroglycan is a transmembrane protein that plays important roles in connecting the extracellular matrix to the cytoskeleton. Acts as a cell adhesion receptor in both muscle and non-muscle tissues. Receptor for both DMD and UTRN and, through these interactions, scaffolds axin to the cytoskeleton. Also functions in cell adhesion-mediated signaling and implicated in cell polarity.</p>
组织特异性	Expressed in a variety of fetal and adult tissues. In epidermal tissue, located to the basement membrane. Also expressed in keratinocytes and fibroblasts.
疾病相关	<p>Defects in DAG1 are the cause of muscular dystrophy-dystroglycanopathy limb-girdle type C7 (MDDGC7) [MIM:613818]. An autosomal recessive muscular dystrophy showing onset in early childhood, and associated with mental retardation without structural brain anomalies.</p> <p>Note=MDDGC7 is caused by DAG1 mutations that interfere with normal post-translational processing, resulting in defective DAG1 glycosylation and impaired interactions with extracellular-matrix components. Other muscular dystrophy-dystroglycanopathies are caused by defects in enzymes involved in protein O-glycosylation.</p>
序列相似性	Contains 1 peptidase S72 domain.
翻译后修饰	O- and N-glycosylated. Alpha-dystroglycan is heavily O-glycosylated comprising of up to two thirds of its mass and the carbohydrate composition differs depending on tissue type. Mucin-type O-glycosylation is important for ligand binding activity. O-mannosylation of alpha-DAG1 is found in high abundance in both brain and muscle where the most abundant glycan is Sia-alpha-2-3-Gal-beta-1-4-Glc-NAc-beta-1-2-Man. In muscle, glycosylation on Thr-379 by a phosphorylated O-mannosyl glycan with the structure 2-(N-acetylamido)-2-deoxygalactosyl-beta-1,3-2-(N-acetylamido)-2-deoxyglucosyl-beta-1,4-6-phosphomannose is mediated by like-

acetylglucosaminyltransferase (LARGE) protein and is required for laminin binding. O-mannosylation is also required for binding lymphocytic choriomeningitis virus, Old World Lassa fever virus, and clade C New World arenaviruses. The O-glycosyl hexose on Thr-367, Thr-369, Thr-372, Thr-381 and Thr-388 is probably mannose. O-glycosylated in the N-terminal region with a core 1 or possibly core 8 glycan. The beta subunit is N-glycosylated.


Autolytic cleavage produces the alpha and beta subunits. In cutaneous cells, as well as in certain pathological conditions, shedding of beta-dystroglycan can occur releasing a peptide of about 30 kDa.

SRC-mediated phosphorylation of the PPXY motif of the beta subunit recruits SH2 domain-containing proteins, but inhibits binding to WWW domain-containing proteins, DMD and UTRN. This phosphorylation also inhibits nuclear entry.

细胞定位

Secreted > extracellular space and Cell membrane. Cytoplasm > cytoskeleton. Nucleus > nucleoplasm. The monomeric form translocates to the nucleus via the action of importins and depends on RAN. Nuclear transport is inhibited by Tyr-892 phosphorylation. In skeletal muscle, this phosphorylated form locates to a vesicular internal membrane compartment. In peripheral nerves, localizes to the Schwann cell membrane. Colocalizes with ERM proteins in Schwann-cell microvilli.

图片



All lanes : Anti-Alpha Dystroglycan antibody [2238] (ab106110) at 1/100 dilution

Lane 1 : Human brain tissue lysate

Lane 2 : SH-SY5Y cell lysate

Lane 3 : Mouse brain tissue lysate

Lane 4 : Mouse heart tissue lysate

Lane 5 : Rat brain tissue lysate

Lane 6 : Rat heart tissue lysate

Lysates/proteins at 20 µg per lane.

Performed under reducing conditions.

Predicted band size: 97 kDa

Observed band size: 75-100 kDa

Lanes 1 - 6: Merged signal (red and green). Green - ab106110 observed at 75-100 kDa. Red - loading control **ab181602** (Rabbit Anti-GAPDH antibody [EPR16891]) observed at 37 kDa.

ab106110 was shown to react with DAG1 in Western blot.

Membranes were blocked in fluorescent western blot (TBS-based) blocking solution before incubation with ab106110 and **ab181602**

(Rabbit Anti-GAPDH antibody [EPR16891]) overnight at 4 °C at a 1 in 100 dilution and a 1 in 20000 dilution respectively. Blots were incubated with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed (**ab216773**) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed (**ab216776**) secondary antibodies at 1 in 20000 dilution for 1 h at room temperature before imaging.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors