

Anti-ALPK2 antibody ab111909

5 References [1 图像](#)

概述

产品名称	Anti-ALPK2抗体
描述	兔多克隆抗体to ALPK2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	A peptide from C terminal residues of Human ALPK2 protein (NM_052947.3).
阳性对照	Human fetal liver lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Lyophilized:Reconstitute in 200ul Sterile Distilled Water.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 1% BSA, 98.98% PBS</p>
纯度	Immunogen affinity purified
纯化说明	ab111909 is purified by a peptide affinity column.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab111909于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 237 kDa (predicted molecular weight: 237 kDa).

靶标

功能	Kinases that recognize phosphorylation sites in which the surrounding peptides have an alpha-helical conformation.
序列相似性	Belongs to the protein kinase superfamily. Alpha-type protein kinase family. ALPK subfamily. Contains 1 alpha-type protein kinase domain. Contains 2 Ig-like (immunoglobulin-like) domains.

图片



Anti-ALPK2 antibody (ab111909) at 1/500 dilution + Human fetal liver lysate

Predicted band size: 237 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors