

APC Anti-Ly76 antibody [TER-119] ab93582

1 图像

概述

产品名称	APC Anti-Ly76抗体[TER-119]
描述	APC大鼠单克隆抗体[TER-119] to Ly76
宿主	Rat
偶联物	APC. Ex: 645nm, Em: 660nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue/ cell preparation: Mouse fetal liver cells.
阳性对照	Mouse erythroid cells from early proerythroblast to mature erythrocyte stages. Mouse splenocytes and bone marrow cell suspensions.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 0.87% Sodium chloride, 0.14% Monobasic dihydrogen sodium phosphate, 0.1% Gelatin</p>
纯度	Affinity purified
克隆	单克隆
克隆编号	TER-119
同种型	IgG2b
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab93582于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration.

靶标

相关性 Ly76 is a glycoprotein A-associated protein. It is an erythroid-specific antigen expressed on early proerythroblasts to mature erythrocytes, but not BFU-E or CFU-E.

细胞定位 Plasma membrane

图片



Flow cytometry of C57Bl/6 bone marrow cells labeling Ly76 with ab93582 at 0.5 ug (right panel) or Rat IgG2b isotype control at 0.5 ug (left panel).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors