

# APC Mouse IgG1 [1F8] - Isotype Control ab91358

## 2 References

### 概述

产品名称	APC小鼠IgG1 [1F8] -同型对照
偶联物	APC. Ex: 645nm, Em: 660nm
特异性	ab91358 was produced against a synthetic hapten, which is normally not present in Humans or animals. Therefore, this mouse IgG1 does not react with any antigen of Human or animal cells.
经测试应用	<b>适用于:</b> Flow Cyt
免疫原	Vegetal protein
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.20 Preservative: 0.1% Sodium azide Constituents: PBS, 1% BSA
纯度	Protein G purified
克隆	单克隆
克隆编号	1F8
同种型	IgG1

### 应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab91358于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		

#### 应用说明

Flow Cyt: Use 20µl for 10<sup>6</sup> cells.  
ab91358 is a negative isotype control

Not yet tested in other applications.  
Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors