

APC Anti-NCAM1 antibody [LT56] ab205500

1 References [1 图像](#)

概述

产品名称	APC Anti-NCAM1抗体[LT56]
描述	APC小鼠单克隆抗体[LT56] to NCAM1
宿主	Mouse
偶联物	APC. Ex: 645nm, Em: 660nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human NCAM. Cell line KG1a.
阳性对照	Human blood cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	<p>pH: 7.4</p> <p>Preservative: 0.0975% Sodium azide</p> <p>Constituents: 99% PBS, 0.2% BSA</p>
纯度	High grade protease free BSA.
克隆	Size exclusion
克隆编号	单克隆
同种型	LT56
	IgG2a

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab205500于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Use 10 µl for 100 µl of Human whole blood or 10 ⁶ cells in suspension.

靶标

功能

This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.

序列相似性


Contains 2 fibronectin type-III domains.

Contains 5 Ig-like C2-type (immunoglobulin-like) domains.

细胞定位

Secreted and Cell membrane.

图片



Flow Cytometry analysis of human peripheral blood cells labeling NULL with Anti-NCAM antibody [LT56] (Allophycocyanin) (ab205500).

Flow Cytometry - APC Anti-NCAM1 antibody [LT56]
(ab205500)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors