

Anti-ASF1 A+B antibody ab53608

[1 Abreviews](#) [1 References](#) [1 图像](#)

概述	
产品名称	Anti-ASF1 A+B抗体
描述	兔多克隆抗体to ASF1 A+B
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Drosophila melanogaster
免疫原	Recombinant full length protein (Fruit fly (Drosophila melanogaster))
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: 50% Glycerol, Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	Abpromise™ 承诺保证使用ab53608于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	


应用	Ab评论	说明
WB		1/2000. Predicted molecular weight: 24 kDa.

靶标

相关性 Asf1 (anti silencing function 1) belongs to the H3/H4 family of histone chaperones. Its function has not been fully elucidated but it is believed to be involved in ensuring a constant supply of histones at sites of nucleosome assembly. Asf1 interacts with the DNA damage checkpoint protein Rad53 and so may have a role in cell cycle regulation and DNA responses. ASF1 A and B (the two human subtypes) are the only known substrate for Tousled like kinase (TLK).

细胞定位 Nuclear

图片



All lanes : Anti-ASF1 A+B antibody (ab53608) at 1/2000 dilution

All lanes : Drosophila embryo lysate

Predicted band size: 24 kDa

Observed band size: 30 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors