

### Anti-AUTS2 antibody ab96326

★★★★★ [1 Abreviews](#) [4 References](#) [2 图像](#)

#### 概述

产品名称	Anti-AUTS2抗体
描述	兔多克隆抗体to AUTS2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 524 and 770 of AUTS2
阳性对照	Raji whole cell lysate ( <a href="#">ab30124</a> ) can be used as a positive control in WB. 293T and NIH-3T3 cell lysate
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.00 Preservative: 0.025% Proclin 300 Constituents: 79% PBS, 20% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab96326 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	1/500 - 1/3000. Predicted molecular weight: 139 kDa.

## 靶标


### 组织特异性

Strongly expressed in brain, skeletal muscle and kidney. Also expressed in placenta, lung and leukocytes.

### 序列相似性


Belongs to the AUTS2 family.

## 图片



Anti-AUTS2 antibody (ab96326) at 1/1000 dilution + 293T cell lysate at 30  $\mu$ g

**Predicted band size:** 139 kDa



Anti-AUTS2 antibody (ab96326) at 1/1000 dilution + NIH3T3 cell lysate at 30  $\mu$ g

**Predicted band size:** 139 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors