

Anti-Benzoylecgonine antibody [3P1A6] ab116601

概述

产品名称	Anti-Benzoylecgonine抗体[3P1A6]
描述	小鼠单克隆抗体[3P1A6] to Benzoylecgonine
宿主	Mouse
经测试应用	适用于: IP, ELISA, RIA
种属反应性	与反应: Species independent
免疫原	Benzoylecgonine conjugated to KLH.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
存储溶液	Preservative: 0.05% Sodium azide Constituents: 0.019% Potassium chloride, 0.14% Potassium phosphate, 0.79% Sodium chloride, 0.1% BSA
纯度	Protein A purified
克隆	单克隆
克隆编号	3P1A6
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab116601于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
RIA		Use at an assay dependent dilution.

靶标

相关性 Benzoylcegonine is the predominant cocaine metabolite. This compound is not biologically active but is used as an indicator of cocaine use.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors