

Biotin Anti-CD90 / Thy1 antibody [5E10] ab134361

[1 References](#) [1 图像](#)

概述

产品名称	生物素Anti-CD90 / Thy1抗体[5E10]
描述	生物素小鼠单克隆抗体[5E10] to CD90 / Thy1
宿主	Mouse
偶联物	Biotin
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human CD90/ Thy1. HEL Human erythroleukemia cell line
阳性对照	Human erythroleukemia (HEL) cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at +4°C. Do Not Freeze. Store In the Dark.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.1% Gelatin, 99% PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	5E10
同种型	IgG1
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用 ab134361 于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 1 µg for 10 ⁵⁻⁸ cells. (In a final volume of 100 µl.) ab134361 Mouse monoclonal IgG1 is suitable for use as an

靶标

功能	May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.
序列相似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
细胞定位	Cell membrane.

图片



Flow cytometry analysis of Human erythroleukemia (HEL) cells with 0.5 µg of Mouse IgG1 kappa Isotype Control Biotin (open histogram) or 0.5 µg ab134361 (filled histogram) followed by Streptavidin PE. Total viable cells were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors