

### Anti-C5a-R antibody [10/92] ab117579

**4 References**   [1 图像](#)

#### 概述

产品名称	Anti-C5a-R抗体[10/92]
描述	大鼠单克隆抗体[10/92] to C5a-R
宿主	Rat
经测试应用	适用于: IHC-P
种属反应性	与反应: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse C5a-R. RBL-2H3 transfected cells
阳性对照	Mouse Kidney tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: 99% PBS</p>
纯度	Protein G purified
克隆	单克隆
克隆编号	10/92
同种型	IgG2a

#### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab117579于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 10 µg/ml.

## 靶标

功能	Receptor for the chemotactic and inflammatory peptide anaphylatoxin C5a. This receptor stimulates chemotaxis, granule enzyme release and superoxide anion production.
序列相似性	Belongs to the G-protein coupled receptor 1 family.
翻译后修饰	Sulfation plays a critical role in the association of the receptor with C5a, but no significant role in the ability of the receptor to transduce a signal and mobilize calcium in response to a small peptide agonist.
细胞定位	Cell membrane.

## 图片



Immunohistochemical analysis of C5a-R expression in Formalin-Fixed, Paraffin-Embedded Mouse Kidney tissue using 10 µg/ml ab117579.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C5a-R antibody [10/92] (ab117579)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors