

### Anti-CA19-9 antibody [121SLE] ab3982

★★★★★ [1 Abreviews](#) [16 References](#) [1 图像](#)

#### 概述

产品名称	Anti-CA19-9抗体[121SLE]
描述	小鼠单克隆抗体[121SLE] to CA19-9
宿主	Mouse
特异性	The antibody recognizes Sialyl Lewis a containing glycolipids.
经测试应用	<b>适用于:</b> IHC-P, ELISA, Flow Cyt
种属反应性	<b>与反应:</b> Human
免疫原	Tissue, cells or virus corresponding to CA19-9. Immunoprecipitate obtained after immunodiffusion of MAb 19-9 and mucins isolated from an ovarian cyst of a 0Le(a+b-) patient.
阳性对照	IHC-P: Colorectal adenocarcinoma.
常规说明	<p>This product was changed from ascites to tissue culture supernatant. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze.
存储溶液	<p>pH: 8.0</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Proprietary Purification
纯化说明	Purified from TCS. Purified by precipitation and chromatography. Purity >95% by SDS-PAGE.
克隆	单克隆

克隆编号121SLE

同种型IgM

应用


The Abpromise guarantee      **Abpromise™**承诺保证使用ab3982于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/100.
ELISA		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. PubMed: 22490524 <b>ab91545</b> - Mouse monoclonal IgM, is suitable for use as an isotype control with this antibody.

靶标

图片



ab3982 staining CA19-9 in human colorectal adenocarcinoma by IHC-P (Formalin/Paraformaldehyde fixed paraffin embedded tissue sections).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors