

Anti-CAPON/NOS1AP antibody [EPR14519(B)] ab190686

重组 RabMAb

3 References **6 图像**

概述

产品名称	Anti-CAPON/NOS1AP抗体[EPR14519(B)]
描述	兔单克隆抗体[EPR14519(B)] to CAPON/NOS1AP
宿主	Rabbit
经测试应用	适用于: Flow Cyt (Intra), ICC/IF, WB
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	Human fetal brain, Human cerebellum, Mouse brain and Rat brain lysates: MCF-7, HeLa, C6, Raw 264.7, PC-12 and NIH/3T3 cell lysates: SH-SY5Y and HeLa cells.
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR14519(B)
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab190686于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Flow Cyt (Intra)		1/130. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
ICC/IF		1/50.
WB		1/1000 - 1/10000. Detects a band of approximately 56 kDa (predicted molecular weight: 56 kDa).

靶标

功能 Adapter protein involved in neuronal nitric-oxide (NO) synthesis regulation via its association with nNOS/NOS1. The complex formed with NOS1 and synapsins is necessary for specific NO and synapsin functions at a presynaptic level. Mediates an indirect interaction between NOS1 and RASD1 leading to enhance the ability of NOS1 to activate RASD1. Competes with DLG4 for interaction with NOS1, possibly affecting NOS1 activity by regulating the interaction between NOS1 and DLG4.

序列相似性 Contains 1 PID domain.

图片




All lanes : Anti-CAPON/NOS1AP antibody [EPR14519(B)] (ab190686) at 1/1000 dilution

Lane 1 : Human fetal brain lysates
Lane 2 : Human cerebellum brain lysates
Lane 3 : MCF-7 cell lysates
Lane 4 : HeLa cell lysates

Lysates/proteins at 20 µg per lane.

Secondary
All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 56 kDa




Western blot - Anti-CAPON/NOS1AP antibody [EPR14519(B)] (ab190686)

Anti-CAPON/NOS1AP antibody [EPR14519(B)] (ab190686) at 1/10000 dilution + Human fetal brain lysates at 20 µg

Secondary

Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 56 kDa



Western blot - Anti-CAPON/NOS1AP antibody [EPR14519(B)] (ab190686)

All lanes : Anti-CAPON/NOS1AP antibody [EPR14519(B)] (ab190686) at 1/1000 dilution

Lane 1 : Mouse brain lysates

Lane 2 : Rat brain lysates

Lane 3 : C6 cell lysates

Lane 4 : Raw 264.7 cell lysates

Lane 5 : PC-12 cell lysates


Lane 6 : NIH/3T3 cell lysates

Lysates/proteins at 10 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 56 kDa



Immunocytochemistry/ Immunofluorescence - Anti-CAPON/NOS1AP antibody [EPR14519(B)] (ab190686)

Immunofluorescence analysis of 4% paraformaldehyde-fixed SH-SY5Y cells labeling CAPON/NOS1AP with ab190686 at 1/50 dilution. Goat anti-rabbit IgG (Alexa Fluor® 555) at 1/200 dilution was used as the secondary antibody (red). Slide is counter stained with Dapi (blue).



Flow Cytometry (Intracellular) - Anti-CAPON/NOS1AP antibody [EPR14519(B)] (ab190686)

Intracellular flow cytometric analysis of 2% paraformaldehyde fixed HeLa cells labeling CAPON/NOS1AP with ab190686 at 1/130 dilution (red) compared with a rabbit monoclonal IgG control (green). Goat anti rabbit IgG (FITC) at 1/150 dilution was used as the secondary antibody.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-CAPON/NOS1AP antibody [EPR14519(B)]
(ab190686)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors