

Anti-CD105 antibody [OT18A1] ab156756

★★★★★ [2 Abreviews](#) [15 References](#) [4 图像](#)

概述

产品名称	Anti-CD105抗体[OT18A1]
描述	小鼠单克隆抗体[OT18A1] to CD105
宿主	Mouse
经测试应用	适用于: WB, ICC/IF, Flow Cyt (Intra)
种属反应性	与反应: Mouse, Rat, Dog, Human, African green monkey
免疫原	Recombinant full length protein corresponding to Human CD105.
阳性对照	WB: HepG2, HeLa, SVT2, A549, COS7, Jurkat, MDCK, PC12, and MCF7 cell lysates. IF: HeLa cells. Flow Cyt (Intra): U937 cells.
常规说明	<p>The clone number has been updated from 8A1 to OT18A1, both clone numbers name the same clone.</p> <p>This product was changed from ascites to tissue culture supernatant on 29th May 2018. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 1% BSA, 50% Glycerol
纯度	Immunogen affinity purified
纯化说明	ab156756 is purified from TCS by affinity chromatography.

克隆	单克隆
克隆编号	OT18A1
同种型	IgG2b

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab156756于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/200 - 1/4000. Predicted molecular weight: 70 kDa.
ICC/IF		1/100.
Flow Cyt (Intra)		1/1000. ab170192 - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.

靶标


功能	Major glycoprotein of vascular endothelium. May play a critical role in the binding of endothelial cells to integrins and/or other RGD receptors.
组织特异性	Endoglin is restricted to endothelial cells in all tissues except bone marrow.
疾病相关	Defects in ENG are the cause of hereditary hemorrhagic telangiectasia type 1 (HHT1) [MIM:187300, 108010]; also known as Osler-Rendu-Weber syndrome 1 (ORW1). HHT1 is an autosomal dominant multisystemic vascular dysplasia, characterized by recurrent epistaxis, muco-cutaneous telangiectases, gastro-intestinal hemorrhage, and pulmonary (PAVM), cerebral (CAVM) and hepatic arteriovenous malformations; all secondary manifestations of the underlying vascular dysplasia. Although the first symptom of HHT1 in children is generally nose bleed, there is an important clinical heterogeneity.
细胞定位	Membrane.

图片




Flow Cytometry (Intracellular) - Anti-CD105 antibody [OT18A1] (ab156756)

Overlay histogram showing U937 cells stained with ab156756 (red line). The cells were fixed with 4% paraformaldehyde (10 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab156756, 1/1000 dilution) for 30 min at 22°C. The secondary antibody used was Alexa Fluor® 488 goat anti-mouse IgG (H&L) (**ab150113**) at 1/2000 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG2b [PLPV219] (**ab91366**, 1µg/1x10⁶ cells) used under the same conditions. Unlabelled sample (blue line) was also used as a control. Acquisition of >5,000 events were collected using a 20mW Argon ion laser (488nm) and 525/30 bandpass filter. This antibody gave a positive signal in U937 cells fixed with 80% methanol (5 min)/permeabilized with 0.1% PBS-Tween for 20 min used under the same conditions.



Immunocytochemistry/ Immunofluorescence - Anti-CD105 antibody [OT18A1] (ab156756)

Immunofluorescent analysis of HeLa cells labeling CHN1 with ab156756 at 1/100.



Western blot - Anti-CD105 antibody [OT18A1] (ab156756)

All lanes : Anti-CD105 antibody [OT18A1] (ab156756) at 1/200 dilution


Lane 1 : HEK293T cells transfected with pCMV6-ENTRY control

Lane 2 : HEK293T cells transfected with pCMV6-ENTRY CD105

Lysates/proteins at 5 µg per lane.

Predicted band size: 70 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC211854)



All lanes : Anti-CD105 antibody [OT18A1] (ab156756) at 1/200 dilution

Lane 1 : HepG2 cell lysate

Lane 2 : HeLa cell lysate

Lane 3 : SVT2 cell lysate

Lane 4 : A549 cell lysate

Lane 5 : COS7 cell lysate

Lane 6 : Jurkat cell lysate

Lane 7 : MDCK cell lysate

Lane 8 : PC12 cell lysate

Lane 9 : MCF7 cell lysate

Lysates/proteins at 35 µg per lane.

Predicted band size: 70 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC211696)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors