

Anti-CD3 antibody [RM0027-3B19] ab56313

★★★★★ 2 Abreviews 12 References 2 图像

概述

产品名称	Anti-CD3抗体[RM0027-3B19]
描述	大鼠单克隆抗体[RM0027-3B19] to CD3
宿主	Rat
经测试应用	适用于: IHC-P
种属反应性	与反应: Mouse
免疫原	Mouse CD3 positive cells
阳性对照	Spleen tissue
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	RM0027-3B19
同种型	IgG2b


应用

The Abpromise guarantee Abpromise™承诺保证使用ab56313于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
IHC-P	★★★★★ (2)	1/50 - 1/400.

图片



ab56313, staining Mouse CD3 in Kidney, using Immunohistochemistry, Paraffin embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD3 antibody [RM0027-3B19] (ab56313)



ab56313 at 1/400 dilution staining CD3 in LPS treated mouse spleen sections by immunohistochemistry (Formalin/ PFA fixed paraffin-embedded tissue sections).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD3 antibody [RM0027-3B19] (ab56313)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors