

Anti-CD62E antibody [1.2B6] ab6630


★★★★☆ [3 Abreviews](#) [11 References](#) [1 图像](#)

概述

产品名称	Anti-CD62E抗体[1.2B6]
描述	小鼠单克隆抗体[1.2B6] to CD62E
宿主	Mouse
特异性	This antibody recognises the human CD62E and CD62P cell surface antigens binding to a common epitope on these members of the selectin family.
经测试应用	适用于: Flow Cyt 不适用于: IHC-P
种属反应性	与反应: Human
免疫原	Full length protein corresponding to Human CD62E. Database link: P16581
常规说明	<p>This antibody inhibits neutrophil adhesion to E-selectin and to P-selectin. E-selectin (ELAM-1) is expressed on the surface of activated endothelial cells.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein A purified
纯化说明	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant.



ab6630 staining CD62E in thrombin activated human peripheral blood platelets by Flow Cytometry analysis.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors