

Anti-CD99 antibody [SPM596] - BSA and Azide free ab212608

1 图像

概述

产品名称	Anti-CD99抗体[SPM596] - BSA and Azide free
描述	小鼠单克隆抗体[SPM596] to CD99 - BSA and Azide free
宿主	Mouse
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human CD99 aa 1-185. Database link: P14209
阳性对照	MOLT-4 cells. Pancreas or Ewing's sarcoma.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2 Constituent: 99% PBS
无载体	是
纯度	Protein G purified
纯化说明	ab212608 was purified from Bioreactor Concentrate.
克隆	单克隆
克隆编号	SPM596
同种型	IgG1
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab212608于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 1 - 2 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. 30 min at Room Temp.

靶标

功能	Involved in T-cell adhesion processes. It is involved in spontaneous rosette formation with erythrocytes.
序列相似性	Belongs to the CD99 family.
翻译后修饰	Extensively O-glycosylated.
细胞定位	Membrane.

图片



Immunohistochemical analysis of formalin-fixed, paraffin-embedded Human Ewing's sarcoma labeling CD99 with ab212608 at 2 µg/ml for 30 minutes at RT.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD99 antibody [SPM596]
- BSA and Azide free (ab212608)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors