

### Anti-CDYL antibody ab5188

★★★★★ [2 Abreviews](#) [13 References](#) [2 图像](#)

#### 概述

产品名称	Anti-CDYL 抗体
描述	兔多克隆抗体 to CDYL
宿主	Rabbit
经测试应用	适用于: IHC-P, ChIP, WB, ChIP-sequencing
种属反应性	与反应: Mouse, Rat, Human, African green monkey
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	This antibody gave a positive signal in COS7 cells transfected with Human CDYL (Whole Cell Lysate)
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: PBS</p> <p>Batches of this product that have a concentration &lt; 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.</p>
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee      **Abpromise™**承诺保证使用ab5188于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
ChIP	★★★★★ (1)	Use at an assay dependent concentration. PubMed: 22009739
WB	★★★★★ (1)	1/1000. Detects a band of approximately 80 kDa (predicted molecular weight: 66 kDa).
ChIP-sequencing		Use at an assay dependent concentration. PubMed: 22009739

靶标

功能	Acts as repressor of transcription (By similarity). Has histone acetyltransferase activity, with a preference for histone H4.
组织特异性	Ubiquitous. Expressed at moderate levels in all tissues examined.
序列相似性	Contains 1 chromo domain.
细胞定位	Nucleus.

图片



**All lanes :** Anti-CDYL antibody (ab5188) at 1 µg/ml

**Lane 1 :** COS-7 (African green monkey kidney fibroblast-like cell line) whole cell lysate - Untransfected

**Lane 2 :** COS7 CDYL (Human) Transfected Whole Cell Lysate

Lysates/proteins at 10 µg per lane.

**Secondary**

**All lanes :** Goat polyclonal to Rabbit IgG - H&L - Pre-Adsorbed (HRP) at 1/3000 dilution

Performed under reducing conditions.


**Predicted band size:** 66 kDa

**Observed band size:** 80 kDa

**Additional bands at:** 60 kDa, 65 kDa. We are unsure as to the

identity of these extra bands.

**Exposure time:** 12 minutes



**All lanes :** Anti-CDYL antibody (ab5188) at 1/1000 dilution

**Lane 1 :** COS-7 (African green monkey kidney fibroblast-like cell line) whole cell lysate

**Lane 2 :** COS7 cells transfected with HA-CDYL-FLAG plasmid

**Predicted band size:** 66 kDa

**Observed band size:** 83 kDa

Rabbit polyclonal to CDYL (ab5188) at 1/1000.

Left lane (-) COS7 (monkey) cells.

Right lane (+) COS7 cells transfected with HA-CDYL-FLAG plasmid.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors