

### Anti-CGA/hCG alpha antibody [LHa/756] ab219296

#### 1 图像

#### 概述

产品名称	Anti-CGA/hCG alpha抗体[LHa/756]
描述	小鼠单克隆抗体[LHa/756] to CGA/hCG alpha
宿主	Mouse
经测试应用	适用于: IHC-P
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human CGA/hCG alpha. Database link: <a href="#">P01215</a>
阳性对照	IHC-P: Human pituitary tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2 Preservative: 0.05% Sodium azide Constituents: 99% PBS, 0.05% BSA
纯度	Protein A/G purified
克隆	单克隆
克隆编号	LHa/756
同种型	IgG1
轻链类型	kappa



## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors