

# Anti-Chloramphenicol Acetyltransferase antibody [CAT-1] ab5410

## 1 References

### 概述

产品名称	Anti-Chloramphenicol Acetyltransferase抗体[CAT-1]
描述	小鼠单克隆抗体[CAT-1] to Chloramphenicol Acetyltransferase
宿主	Mouse
经测试应用	适用于: WB, ICC
种属反应性	与反应: Escherichia coli
免疫原	Recombinant full length protein.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	CAT-1
骨髓瘤	Sp2/0-Ag14
同种型	IgG1

### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab5410于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use a concentration of 2 µg/ml.
ICC		Use at an assay dependent dilution.

## 靶标

---

**相关性** CAT is used as a tag in many systems and is used in many transcription assays. The enzyme is an effector of chloramphenicol resistance in bacteria.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors