

### Anti-Collagen III antibody [EPR17673] ab184993

**重组** RabMAb

**47 References**   **16 图像**

#### 概述

产品名称	Anti-Collagen III抗体[EPR17673]
描述	兔单克隆抗体[EPR17673] to Collagen III
宿主	Rabbit
经测试应用	<b>适用于:</b> Flow Cyt (Intra), ICC/IF, IP, WB
种属反应性	<b>与反应:</b> Mouse, Rat, Human
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: Human fetal liver, fetal kidney, skin, skeletal muscle and heart lysates; HeLa, MCF7, HaCaT and A431 whole cell lysates. Mouse skeletal muscle tissue lysates, rat skeletal muscle tissue lysates, mouse heart tissue lysate, C6, RAW264.7, PC-12 and NIH/3T3 whole cell lysates. C6, RAW264.7 and PC-12 whole cell fresh lysate. ICC/IF: HeLa, MCF7, PC-12 and RAW 264.7 cells Flow Cyt (intra): HeLa, PC-12 and RAW 264.7 cells IP: HeLa whole cell lysate.
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> <li>- High batch-to-batch consistency and reproducibility</li> <li>- Improved sensitivity and specificity</li> <li>- Long-term security of supply</li> <li>- Animal-free production</li> </ul> <p>For more information <a href="#">see here</a>.</p> <p>Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a>.</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 0.05% BSA, 40% Glycerol, PBS</p>
纯度	Protein A purified
克隆	单克隆

克隆编号	EPR17673
同种型	IgG

## 应用


**The Abpromise guarantee** Abpromise™承诺保证使用ab184993于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Flow Cyt (Intra)		1/50.
ICC/IF		1/100.
IP		1/30.
WB		1/1000. Detects a band of approximately 150 kDa (predicted molecular weight: 139 kDa).

## 靶标

功能	Collagen type III occurs in most soft connective tissues along with type I collagen.
疾病相关	<p>Defects in COL3A1 are a cause of Ehlers-Danlos syndrome type 3 (EDS3) [MIM:130020]; also known as benign hypermobility syndrome. EDS is a connective tissue disorder characterized by hyperextensible skin, atrophic cutaneous scars due to tissue fragility and joint hyperlaxity. EDS3 is a form of Ehlers-Danlos syndrome characterized by marked joint hyperextensibility without skeletal deformity.</p> <p>Defects in COL3A1 are the cause of Ehlers-Danlos syndrome type 4 (EDS4) [MIM:130050]. EDS is a connective tissue disorder characterized by hyperextensible skin, atrophic cutaneous scars due to tissue fragility and joint hyperlaxity. EDS4 is the most severe form of the disease. It is characterized by the joint and dermal manifestations as in other forms of the syndrome, characteristic facial features (acrogeria) in most patients, and by proneness to spontaneous rupture of bowel and large arteries. The vascular complications may affect all anatomical areas.</p> <p>Defects in COL3A1 are a cause of susceptibility to aortic aneurysm abdominal (AAA) [MIM:100070]. AAA is a common multifactorial disorder characterized by permanent dilation of the abdominal aorta, usually due to degenerative changes in the aortic wall. Histologically, AAA is characterized by signs of chronic inflammation, destructive remodeling of the extracellular matrix, and depletion of vascular smooth muscle cells.</p>
序列相似性	<p>Belongs to the fibrillar collagen family.</p> <p>Contains 1 fibrillar collagen NC1 domain.</p> <p>Contains 1 VWFC domain.</p>
翻译后修饰	<p>Proline residues at the third position of the tripeptide repeating unit (G-X-Y) are hydroxylated in some or all of the chains.</p> <p>O-linked glycan consists of a Glc-Gal disaccharide bound to the oxygen atom of a post-translationally added hydroxyl group.</p>
细胞定位	Secreted > extracellular space > extracellular matrix.



Immunocytochemistry/ Immunofluorescence - Anti-Collagen III antibody [EPR17673] (ab184993)



Western blot - Anti-Collagen III antibody [EPR17673] (ab184993)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized MCF7 (human breast adenocarcinoma cell line) cells labeling Collagen III with ab184993 at 1/100 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) ([ab150077](#)) secondary antibody at 1/1000 dilution (green). Confocal image showing cytoplasmic staining on MCF7 cells.

The nuclear counter stain is DAPI (blue). Tubulin is detected with Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) ([ab195889](#)) (red) at 1/200 dilution.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) ([ab150077](#)) secondary antibody at 1/1000 dilution.

**All lanes :** Anti-Collagen III antibody [EPR17673] (ab184993) at 1/1000 dilution

- Lane 1 :** Human fetal liver lysate at 20 µg
- Lane 2 :** Human fetal kidney lysate at 20 µg
- Lane 3 :** HeLa (human epithelial cell line from cervix adenocarcinoma) whole cell lysate at 20 µg
- Lane 4 :** MCF7 (human breast adenocarcinoma cell line) whole cell lysate at 20 µg
- Lane 5 :** HaCaT (human keratinocyte cell line) whole cell lysate at 10 µg

**Secondary**

**All lanes :** Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution

Developed using the ECL technique.


**Predicted band size:** 139 kDa  
**Observed band size:** 150 kDa

**Exposure time :** Lane 1: 3 minutes; Lane 2: 15 seconds; Lane 3: 3 seconds; Lanes 4-5: 30 seconds.

**Blocking/Dilution buffer:** 5% NFDM/TBST.

The molecular weight observed is consistent with what has been

described in the literature (PMID: 26017148).



Immunoprecipitation - Anti-Collagen III antibody  
[EPR17673] (ab184993)

Collagen III was immunoprecipitated from 0.35 mg of HeLa (human epithelial cell line from cervix adenocarcinoma) whole cell lysate with ab184993 at 1/30 dilution. Western blot was performed from the immunoprecipitate using ab184993 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) ([ab131366](#)), was used for detection at 1/5000 dilution.


**Lane 1:** HeLa whole cell lysate 10 µg (Input).

**Lane 2:** ab184993 IP in HeLa whole cell lysate .

**Lane 3:** Rabbit monoclonal IgG ([ab172730](#)) instead of ab184993 in HeLa whole cell lysate .

Blocking and dilution buffer and concentration: 5% NFDM/TBST.

Exposure time: 5.5 seconds.




Immunocytochemistry/ Immunofluorescence - Anti-Collagen III antibody [EPR17673] (ab184993)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized HeLa (human epithelial cell line from cervix adenocarcinoma) cells labeling Collagen III with ab184993 at 1/100 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) ([ab150077](#)) secondary antibody at 1/1000 dilution (green). Confocal image showing cytoplasmic staining on HeLa cells.

The nuclear counter stain is DAPI (blue). Tubulin is detected with Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) ([ab195889](#)) (red) at 1/200 dilution.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) ([ab150077](#)) secondary antibody at 1/1000 dilution.




Immunocytochemistry/ Immunofluorescence - Anti-Collagen III antibody [EPR17673] (ab184993)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized PC-12 (rat adrenal gland pheochromocytoma small irregularly shaped cells) cell line labeling Collagen III with ab184993 at 1/50 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) preadsorbed secondary antibody (**ab150081**) at 1/1000 dilution (green). Confocal image showing cytoplasmic staining in PC-12 cell line.

The nuclear counter stain is DAPI (blue). Tubulin is detected with Anti-alpha Tubulin mouse monoclonal antibody - Microtubule Marker (Alexa Fluor® 594) (**ab195889**) (red) at 1/200 dilution.

Secondary antibody only control: Used PBS instead of primary antibody, secondary is Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) preadsorbed secondary antibody (**ab150081**) at 1/1000 dilution.




Immunocytochemistry/ Immunofluorescence - Anti-Collagen III antibody [EPR17673] (ab184993)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized RAW 264.7 (mouse abelson murine leukemia virus-induced tumor macrophage) cell line labeling Collagen III with ab184993 at 1/50 dilution followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) preadsorbed secondary antibody (**ab150081**) at 1/1000 dilution (green). Confocal image showing cytoplasmic staining in RAW 264.7 cell line.

The nuclear counter stain is DAPI (blue). Tubulin is detected with Anti-alpha Tubulin mouse monoclonal antibody - Microtubule Marker (Alexa Fluor® 594) (**ab195889**) (red) at 1/200 dilution.

Secondary antibody only control: Used PBS instead of primary antibody, secondary is Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) preadsorbed secondary antibody (**ab150081**) at 1/1000 dilution.



Immunoprecipitation - Anti-Collagen III antibody [EPR17673] (ab184993)

Collagen III was immunoprecipitated from 0.35 mg of RAW 264.7 (mouse Abelson murine leukemia virus-induced tumor macrophage), whole cell lysate with ab184993 at 1/30 dilution. Western blot was performed from the immunoprecipitate using ab184993 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) ([ab131366](#)), was used for detection at 1/5000 dilution.

**Lane 1:** RAW 264.7 whole cell lysate 10 µg (Input).


**Lane 2:** ab184993 IP in RAW 264.7 whole cell lysate.

**Lane 3:** Rabbit monoclonal IgG ([ab172730](#)) instead of ab184993 in RAW 264.7 whole cell lysate.

Blocking and dilution buffer and concentration: 5% NFDM/TBST.

Exposure time: 3.25 seconds.

Lysate was freshly made and used immediately to minimize protein degradation.



Western blot - Anti-Collagen III antibody [EPR17673] (ab184993)

**All lanes :** Anti-Collagen III antibody [EPR17673] (ab184993) at 1/1000 dilution

**Lane 1 :** Human skin lysate

**Lane 2 :** Human skeletal muscle lysate

**Lane 3 :** Human heart lysate

**Lane 4 :** A431 (human epidermoid carcinoma cell line) whole cell lysate

Lysates/proteins at 10 µg per lane.

### Secondary

**Lane 1 :** Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution

**Lanes 2-4 :** Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Developed using the ECL technique.


**Predicted band size:** 139 kDa

**Observed band size:** 150 kDa

**Exposure time :** Lanes 1-3: 30 seconds; Lane 4: 5 seconds.

**Blocking/Dilution buffer:** 5% NFDM/TBST.

The molecular weight observed is consistent with what has been described in the literature (PMID: 26017148).



**All lanes :** Anti-Collagen III antibody [EPR17673] (ab184993) at 1/1000 dilution

**Lane 1 :** C6 (rat glial tumor glial cell), whole cell lysate

**Lane 2 :** RAW264.7 (mouse Abelson murine leukemia virus-induced tumor macrophage), whole cell lysate

**Lane 3 :** PC-12 (rat adrenal gland pheochromocytoma), whole cell lysate

**Lane 4 :** NIH/3T3 (mouse embryonic fibroblast), whole cell lysate

**Lane 5 :** C6 (rat glial tumor glial cell), whole cell fresh lysate

**Lane 6 :** RAW264.7 (mouse Abelson murine leukemia virus-induced tumor macrophage), whole cell fresh lysate

**Lane 7 :** PC-12 (rat adrenal gland pheochromocytoma), whole cell fresh lysate

Lysates/proteins at 20 µg per lane.

### Secondary

**All lanes :** Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 139 kDa


**Observed band size:** 150 kDa

**Exposure time :** Lanes 1-4: 92 seconds; Lanes 5-6: 48 seconds; Lane 7: 37 seconds.

**Blocking/Dilution buffer:** 5% NFDM/TBST.

In lane 1-4, the lysates were stored at -80° prior to Western Blotting.

In lane 5-7, the lysates were freshly made and used for Western Blotting immediately to minimize protein degradation.



Western blot - Anti-Collagen III antibody [EPR17673]  
(ab184993)

**All lanes :** Anti-Collagen III antibody [EPR17673] (ab184993) at 1/1000 dilution

**Lane 1 :** Mouse skeletal muscle tissue lysate

**Lane 2 :** Mouse heart tissue lysate

**Lane 3 :** Rat skeletal muscle tissue lysate

Lysates/proteins at 20 µg per lane.

### Secondary

**All lanes :** Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Developed using the ECL technique.

Performed under reducing conditions.


**Predicted band size:** 139 kDa

**Observed band size:** 150 kDa

**Exposure time:** 3 minutes

**Exposure time :** 3 minutes

**Blocking/Dilution buffer:** 5% NFDm/TBST.



Immunoprecipitation - Anti-Collagen III antibody  
[EPR17673] (ab184993)

Collagen III was immunoprecipitated from 0.35 mg of PC-12 (rat adrenal gland pheochromocytoma cell), whole cell lysate with ab184993 at 1/30 dilution. Western blot was performed from the immunoprecipitate using ab184993 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) (**ab131366**), was used for detection at 1/5000 dilution.

**Lane 1:** PC-12 whole cell lysate 10 µg (Input).


**Lane 2:** ab184993 IP in PC-12 whole cell lysate.

**Lane 3:** Rabbit monoclonal IgG (**ab172730**) instead of ab184993 in PC-12 whole cell lysate.

Blocking and dilution buffer and concentration: 5% NFDm/TBST.


Exposure time: 3.25 seconds.

Lysate was freshly made and used immediately to minimize protein degradation.




Flow Cytometry (Intracellular) - Anti-Collagen III antibody [EPR17673] (ab184993)

Intracellular flow cytometric analysis of 4% paraformaldehyde-fixed, 90% methanol permeabilized HeLa (human epithelial cell line from cervix adenocarcinoma) cell line labeling Collagen III with ab184993 at 1/50 dilution (red) compared with a Rabbit IgG, monoclonal [EPR25A] - Isotype Control (**ab172730**) (black) and an unlabeled control (cells without incubation with primary antibody and secondary antibody) (blue). Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) (**ab150077**) at 1/2000 dilution was used as the secondary antibody.




Flow Cytometry (Intracellular) - Anti-Collagen III antibody [EPR17673] (ab184993)

Intracellular flow cytometric analysis of 4% paraformaldehyde-fixed, 90% methanol permeabilized RAW 264.7 (Mouse Abelson murine leukemia virus-induced tumor macrophage) cell line labeling Collagen III with ab184993 at 1/500 dilution (red) compared with a Rabbit IgG, monoclonal [EPR25A] - Isotype Control (**ab172730**) (black) and an unlabeled control (cells without incubation with primary antibody and secondary antibody) (blue). Goat Anti-Rabbit IgG (Alexa Fluor® 488) (**ab150081**) at 1/2000 dilution was used as the secondary antibody.



Flow Cytometry (Intracellular) - Anti-Collagen III antibody [EPR17673] (ab184993)

Intracellular flow cytometric analysis of 4% paraformaldehyde-fixed, 90% methanol permeabilized PC-12 (Rat adrenal gland pheochromocytoma) cell line labeling Collagen III with ab184993 at 1/500 dilution (red) compared with a Rabbit IgG, monoclonal [EPR25A] - Isotype Control (**ab172730**) (black) and an unlabeled control (cells without incubation with primary antibody and secondary antibody) (blue). Goat Anti-Rabbit IgG (Alexa Fluor® 488) (**ab150081**) at 1/2000 dilution was used as the secondary antibody.



Immunoprecipitation - Anti-Collagen III antibody [EPR17673] (ab184993)

Collagen III was immunoprecipitated from 0.35 mg of HeLa (human epithelial cell line from cervix adenocarcinoma) whole cell lysate with ab184993 at 1/30 dilution. Western blot was performed from the immunoprecipitate using ab184993 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) (**ab131366**), was used for detection at 1/10000 dilution.

**Lane 1:** HeLa whole cell lysate 10 µg (Input).


**Lane 2:** ab184993 IP in HeLa whole cell lysate .

**Lane 3:** Rabbit monoclonal IgG (**ab172730**) instead of ab184993 in HeLa whole cell lysate .

Blocking and dilution buffer and concentration: 5% NFDM/TBST.

Exposure time: 10 seconds.

### Why choose a recombinant antibody?



Anti-Collagen III antibody [EPR17673] (ab184993)

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors