

Anti-CXCL9 antibody ab118543

2 图像

概述

产品名称	Anti-CXCL9抗体
描述	兔多克隆抗体to CXCL9
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human
免疫原	Synthetic peptide corresponding to Human CXCL9 aa 1-100 (N terminal). Database link: Q07325

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	Preservatives: 0.025% Thimerosal (merthiolate), 0.01% Sodium azide Constituents: 5% BSA, 0.9% Sodium chloride, 0.1% Dibasic monohydrogen sodium phosphate
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab118543于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use a concentration of 0.5 - 1 µg/ml.
WB		Use a concentration of 0.1 - 0.5 µg/ml. Predicted molecular weight: 14 kDa. The detection limit for CXCL9 is approximately 2.5ng/lane under non-reducing and reducing conditions.

靶标

功能	Cytokine that affects the growth, movement, or activation state of cells that participate in immune and inflammatory response. Chemotactic for activated T-cells. Binds to CXCR3.
序列相似性	Belongs to the intercrine alpha (chemokine CxC) family.
细胞定位	Secreted.


图片



All lanes : Anti-CXCL9 antibody (ab118543)

- Lane 1 :** Recombinant Human CXCL9 Protein at 0.01 µg
- Lane 2 :** Recombinant Human CXCL9 Protein at 0.005 µg
- Lane 3 :** Recombinant Human CXCL9 Protein at 0.0025 µg

Predicted band size: 14 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) staining CXCL9 with Human Lung Cancer Tissue with ab118543

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CXCL9 antibody (ab118543)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors