

Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker ab7754



★★★★★ [8 Abreviews](#) [50 References](#) [9 图像](#)

概述

产品名称	Anti-Cytokeratin 19抗体[A53-B/A2] - Cytoskeleton Marker
描述	小鼠单克隆抗体[A53-B/A2] to Cytokeratin 19 - Cytoskeleton Marker
宿主	Mouse
特异性	Rod domain of cytokeratin peptide 19 (40 kDa) in human tissue.
经测试应用	适用于: ICC/IF, WB, IHC-P, Flow Cyt (Intra)
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human Cytokeratin 19. Human mammary carcinoma cell line MCF-7
表位	Rod domain of cytokeratin peptide 19.
阳性对照	ICC/IF KO: MCF7 cells (MCF7-KRT19 KO used as a negative cell line). HepG2 cells. IHC-P: Human skin. Human liver tissue. WB: HT-29 lysate. Flow Cyt (Intra): MCF7 cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze.
存储溶液	pH: 7.40 Preservative: 0.097% Sodium azide Constituent: PBS
纯度	Protein A purified

克隆	单克隆
克隆编号	A53-B/A2
同种型	IgG2a

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7754于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF	★★★★★ (1)	Use at an assay dependent concentration. Signal can be observed in cells fixed with either methanol or paraformaldehyde.
WB	★★★★★ (3)	Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 44 kDa.
IHC-P	★★★★★ (1)	Use a concentration of 5 - 10 µg/ml. Perform enzymatic antigen retrieval before commencing with IHC staining protocol.
Flow Cyt (Intra)		Use a concentration of 1 - 5 µg/ml. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

靶标


功能	Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.
组织特异性	Expressed in a defined zone of basal keratinocytes in the deep outer root sheath of hair follicles. Also observed in sweat gland and mammary gland ductal and secretory cells, bile ducts, gastrointestinal tract, bladder urothelium, oral epithelia, esophagus, ectocervical epithelium (at protein level). Expressed in epidermal basal cells, in nipple epidermis and a defined region of the hair follicle. Also seen in a subset of vascular wall cells in both the veins and artery of human umbilical cord, and in umbilical cord vascular smooth muscle. Observed in muscle fibers accumulating in the costameres of myoplasm at the sarcolemma in structures that contain dystrophin and spectrin.
序列相似性	Belongs to the intermediate filament family.
发展阶段	Present in hair follicles at all stages of development.
结构域	This keratin differs from all other IF proteins in lacking the C-terminal tail domain.

图片




Flow Cytometry (Intracellular) - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

Separation of MCF-7 cells (red-filled) from human leukocytes (black-dashed) in flow cytometry analysis (intracellular staining) of peripheral whole blood spiked with MCF-7 cells stained using ab7754 (concentration in sample 3 µg/ml, GAM APC).



Immunocytochemistry/ Immunofluorescence - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

ab7754 staining Cytokeratin 19 in wild-type MCF7 cells (top panel) and KRT19 knockout MCF7 cells (bottom panel). The cells were fixed with 100% methanol (5 min), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated with ab7754 at 1/500 dilution and **ab6046** (Rabbit polyclonal to beta Tubulin) at 1/1000 dilution overnight at +4°C, followed by a further incubation at room temperature for 1h with a goat secondary antibody to mouse IgG (Alexa Fluor® 488) (**ab150117**) at 2 µg/ml (shown in green) and a goat secondary antibody to rabbit IgG (Alexa Fluor® 594) (**ab150080**) at 2 µg/ml (shown in red). Nuclear DNA was labelled in blue with DAPI. Image was taken with a confocal microscope (Leica-Microsystems, TCS SP8).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

ab7754 staining Cytokeratin 19 in human skin.

Left panel: with primary antibody at 2 µg/ml. Right panel: isotype control.


Sections were stained using an automated system (DAKO Autostainer Plus), at room temperature: sections were rehydrated and antigen retrieved with the Dako 3 in 1 AR buffers citrate pH6.1. Slides were peroxidase blocked in 3% H₂O₂ in methanol for 10 mins. They were then blocked with Dako Protein block for 10 minutes (containing casein 0.25% in PBS) then incubated with primary antibody for 20 min and detected with Dako envision flex amplification kit for 30 minutes. Colorimetric detection was

completed with Diaminobenzidine for 5 minutes. Slides were counterstained with Haematoxylin and coverslipped under DePeX. Please note that for manual staining we recommend to optimize the primary antibody concentration and incubation time (overnight incubation), and amplification may be required.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

Immunohistochemical analysis of paraffin-embedded human liver tissue stained for Cytokeratin 19 using ab7754 at a 1/100 dilution followed by a GAM IgG-Alexa Fluor[®]488 diluted at 1/200 (green). Cell nuclei stained with PI (1 µg/ml; orange).




Western blot - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

Western blotting analysis of human cytoke... at 2 µg/ml on lysates of HT-29 (Human colorectal adenocarcinoma cell line) cell line and Jurkat (Human T cell leukemia cell line from peripheral blood) cell line (cytokeratin non-expressing cell line; negative control) under non-reducing and reducing conditions.


IRDye800-conjugated anti-mouse IgG1 secondary antibody.

A specific band was detected for cytokeratin 19 at approximately 40 kDa.




Flow Cytometry (Intracellular) - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

Overlay histogram showing MCF7 cells stained with ab7754 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab7754, 0.5µg/1x10⁶ cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG2a [ICIGG2A] ([ab91361](#), 1µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.



Immunocytochemistry/ Immunofluorescence - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)


ICC/IF image of ab7754 stained HepG2 cells. The cells were 100% methanol fixed (5 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody ab7754 at 5µg/ml overnight at +4°C. The secondary antibody (green) was DyLight® 488 goat anti- mouse ([ab96879](#)) IgG (H+L) used at a 1/1000 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

This image is courtesy of an anonymous Abreview

ab7754 staining Cytokeratin 19 in Human liver tissue sections by Immunohistochemistry (IHC-P - formaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 5% milk for 30 minutes at 37°C; antigen retrieval was by heat mediation in citrate buffer. Samples were incubated with primary antibody (1/1000 in antibody diluent) for 1 hour at 37°C. An undiluted HRP-conjugated goat anti-rabbit IgG polyclonal was used as the secondary antibody.



Immunofluorescence analysis of 1 month hepato-differentiated Human dental pulp stem cells, staining Cytokeratin 19 with ab7754 at 1/60 dilution. A FITC-conjugated anti-mouse IgG was used as the secondary antibody.

Immunocytochemistry/ Immunofluorescence - Anti-Cytokeratin 19 antibody [A53-B/A2] - Cytoskeleton Marker (ab7754)

Image from Ferro F et al., PLoS One. 2012;7(7):e41774. Epub 2012 Jul 23. Fig 3.; doi:10.1371/journal.pone.0041774; July 23, 2012, PLoS ONE 7(7): e41774.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors