

Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker ab68459

敲除验证 重组 RabMAb


★★★★★ 1 Abreviews 31 References 17 图像

概述

产品名称	Anti-Cytokeratin 7抗体[EPR1619Y] - Cytoskeleton Marker
描述	兔单克隆抗体[EPR1619Y] to Cytokeratin 7 - Cytoskeleton Marker
宿主	Rabbit
经测试应用	适用于: Flow Cyt (Intra), WB, IP, IHC-P, ICC/IF
种属反应性	与反应: Human
免疫原	Synthetic peptide within Human Cytokeratin 7 aa 1-100 (N terminal). The exact sequence is proprietary.
阳性对照	WB: HeLa, SK-OV-3 and T47D cell lysates. IHC-P: Human ovarian carcinoma and human bladder carcinoma tissues. ICC/IF: HeLa and TD47 cells. Flow Cyt (intra): HeLa cells. IP: HeLa whole cell lysate (ab150035).
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
解离常数 (K _D)	K _D = 2.10 x 10 ⁻¹⁰ M



[Learn more about Kp](#)

存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR1619Y
同种型	IgG

应用

The Abpromise guarantee [Abpromise™](#)承诺保证使用ab68459于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt (Intra)		1/20 - 1/30. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
WB		1/5000 - 1/10000. Detects a band of approximately 51 kDa (predicted molecular weight: 51 kDa).
IP		1/20 - 1/30.
IHC-P		1/1000. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. See IHC antigen retrieval protocols . For unpurified use at 1/100 - 1/250.
ICC/IF	★★★★★ (1)	1/100 - 1/250.

靶标

功能	Blocks interferon-dependent interphase and stimulates DNA synthesis in cells. Involved in the translational regulation of the human papillomavirus type 16 E7 mRNA (HPV16 E7).
组织特异性	Expressed in cultured epidermal, bronchial and mesothelial cells but absent in colon, ectocervix and liver. Observed throughout the glandular cells in the junction between stomach and esophagus but is absent in the esophagus.
序列相似性	Belongs to the intermediate filament family.
翻译后修饰	Arg-20 is dimethylated, probably to asymmetric dimethylarginine.
细胞定位	Cytoplasm.

图片



Western blot - Anti-Cytokeratin 7 antibody
[EPR1619Y] - Cytoskeleton Marker (ab68459)

All lanes : Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459) at 1/5000 dilution

Lane 1 : Wild-type A549 (Human lung carcinoma cell line) whole cell lysate

Lane 2 : KRT7 knockout A549 (Human lung carcinoma cell line) whole cell lysate

Lane 3 : HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate

Lane 4 : MCF7 (Human breast adenocarcinoma cell line) whole cell lysate


Lysates/proteins at 20 µg per lane.

Predicted band size: 51 kDa

Observed band size: 50 kDa

Lanes 1 - 4: Merged signal (red and green). Green - ab68459 observed at 50 kDa. Red - loading control, **ab8245**, observed at 37 kDa.


ab68459 was shown to specifically react with KRT7 (Cytokeratin 7) in wild-type A549 cells as signal was lost in KRT7 knockout cells. Wild-type and KRT7 knockout samples were subjected to SDS-PAGE. Ab68459 and **ab8245** (Mouse anti-GAPDH loading control) were incubated overnight at 4°C at 1/5000 dilution and 1/20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed **ab216773** and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed **ab216776** secondary antibodies at 1/20000 dilution for 1 hour at room temperature before imaging.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human breast carcinoma tissue labelling Cytokeratin 7 with unpurified ab68459.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.




Immunocytochemistry/ Immunofluorescence - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Immunocytochemistry/Immunofluorescence analysis of T47D cells labelling Cytokeratin 7 with purified ab68459 at a dilution of 1/100. Cells were fixed with 4% paraformaldehyde and permeabilized with 0.1% Triton X-100. **ab150077**, an Alexa Fluor® 488-conjugated goat anti-rabbit IgG (1/1000) was used as the secondary antibody. DAPI (blue) was used as the nuclear counterstain. **ab7291**, a mouse anti-tubulin (1/1000) and **ab150120**, an Alexa Fluor® 594-conjugated goat anti-mouse IgG (1/1000) were also used.

Control 1: primary antibody (1/100) and secondary antibody, **ab150120**, an Alexa Fluor® 594-conjugated goat anti-mouse IgG (1/1000).

Control 2: **ab7291** (1/1000) and secondary antibody, **ab150077**, an Alexa Fluor® 488-conjugated goat anti-rabbit IgG (1/1000).

Alexa Fluor® 488 (**ab185048**) and Alexa Fluor® 647 (**ab192077**) conjugated versions are available for this clone.



Western blot - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

All lanes : Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459) at 1/10000 dilution (purified)

Lane 1 : HeLa whole cell lysate

Lane 2 : SK-OV-3 whole cell lysate

Lane 3 : T47D whole cell lysate

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 51 kDa


Observed band size: 51 kDa

Blocking and dilution buffer: 5% NFDM/TBST.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)


Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human bladder carcinoma tissue labelling Cytokeratin 7 with purified ab68459 at a dilution of 1/1000. Heat mediated antigen retrieval was performed using EDTA buffer pH 9. **ab97051**, a HRP-conjugated goat anti-rabbit IgG (H+L) was used as the secondary antibody (1/500). Negative control using PBS instead of primary antibody. Counterstained with hematoxylin.



Flow Cytometry (Intracellular) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Intracellular Flow Cytometry analysis of HeLa cells labelling Cytokeratin 7 with purified ab68459 at a dilution of 1/20 (red). Cells were fixed with 80% methanol. An Alexa Fluor®488-conjugated goat anti-rabbit IgG (1/500) was used as the secondary antibody. Black - Isotype control, rabbit monoclonal IgG. Blue - Unlabelled control, cells without incubation with primary and secondary antibodies.

Alexa Fluor®488 (**ab185048**) and Alexa Fluor®647 (**ab192077**) conjugated versions are available for this clone.



ab68459 (purified) at 1/20 immunoprecipitating Cytokeratin 7 in HeLa whole cell lysate.

Lane 1 (input): HeLa whole cell lysate (10µg)


Lane 2 (+): ab68459 + HeLa whole cell lysate.

Lane 3 (-): Rabbit monoclonal IgG (**ab172730**) instead of ab68459 in HeLa whole cell lysate.

For western blotting, VeriBlot for IP Detection Reagent (HRP) (**ab131366**), was used for detection at 1/1000 dilution.

Blocking buffer and concentration: 5% NFDM/TBST.

Diluting buffer and concentration: 5% NFDM /TBST.



All lanes : Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459) at 1/10000 dilution (unpurified)

Lane 1 : HeLa cell lysate

Lane 2 : T47D cell lysate


Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP-conjugated goat anti-rabbit IgG at 1/2000 dilution

Predicted band size: 51 kDa

Observed band size: 51 kDa



Immunocytochemistry/ Immunofluorescence - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Immunocytochemistry/Immunofluorescence analysis of HeLa cells labelling Cytokeratin 7 with unpurified ab68459 at a dilution of 1/100.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human endometrial carcinoma tissue labelling Cytokeratin 7 with unpurified ab68459.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human ovarian carcinoma tissue labelling Cytokeratin 7 with unpurified ab68459.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human urinary bladder transitional carcinoma tissue labelling Cytokeratin 7 with unpurified ab68459.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

ab68459 showing negative staining in human sarcoma tissue.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

ab68459 showing negative staining in human colonic adenocarcinoma tissue.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Flow Cytometry (Intracellular) - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Overlay histogram showing HeLa cells stained with unpurified ab68459 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (unpurified ab68459, 1/100 dilution) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-rabbit IgG (H+L) ([ab96899](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was rabbit IgG (monoclonal) (1µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.

OL-RD Scanning - Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Equilibrium dissociation constant (K_D)

Learn more about K_D

[Click here to learn more about \$K_D\$](#)

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-Cytokeratin 7 antibody [EPR1619Y] - Cytoskeleton Marker (ab68459)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors