

Anti-Cytokeratin 8+18 antibody [5D3] ab17139

★★★★★ [1 Abreviews](#) [7 References](#) [2 图像](#)

概述

产品名称	Anti-Cytokeratin 8+18抗体[5D3]
描述	小鼠单克隆抗体[5D3] to Cytokeratin 8+18
宿主	Mouse
经测试应用	适用于: WB, ICC/IF, IHC-P, Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue, cells or virus corresponding to Human Cytokeratin 8+18. Cytokeratins from the human breast carcinoma cell line MCF-7
阳性对照	IHC-P: Human Skin tissue. Flow cyt: MCF7 cells.
常规说明	

This antibody stains simple and glandular epithelium. It can be useful for the identification of adenocarcinomas and most squamous cell carcinomas.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.3 Preservative: 0.05% Sodium azide Constituents: Tissue culture supernatant, 1% BSA
纯度	Tissue culture supernatant
Primary antibody说明	This antibody stains simple and glandular epithelium. It can be useful for the identification of adenocarcinomas and most squamous cell carcinomas.

克隆	单克隆
克隆编号	5D3
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用 ab17139 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。


应用	Ab 评论	说明
WB		Use at an assay dependent concentration.
ICC/IF	★★★★★ (1)	Use at an assay dependent concentration. PubMed: 21102523
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
Flow Cyt		1/100. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标

相关性 Function: Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle. Tissue specificity: Observed in muscle fibers accumulating in the costameres of myoplasm at the sarcolemma membrane in structures that contain dystrophin and spectrin. Expressed in gingival mucosa and hard palate of the oral cavity. Disease: Cirrhosis
Similarity: Belongs to the intermediate filament family. PTM: Phosphorylation on serine residues is enhanced during EGF stimulation and mitosis. Ser-74 phosphorylation plays an important role in keratin filament reorganization. O-glycosylated. O-GlcNAcylation at multiple sites increases solubility, and decreases stability by inducing proteasomal degradation. O-glycosylated (O-GlcNAcyated), in a cell cycle-dependent manner.


细胞定位 Cytoplasmic

图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Cytokeratin 8+18 antibody [5D3] (ab17139)

Paraffin embedded human skin tissue stained for Cytokeratin 8/18 using ab 17139 in immunohistochemical analysis.



Flow Cytometry - Anti-Cytokeratin 8+18 antibody [5D3] (ab17139)

Overlay histogram showing MCF7 cells stained with ab17139 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab17139, 1/100 dilution) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG1 [ICIGG1] ([ab91353](#), 2µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed. This antibody gave a positive signal in MCF7 cells fixed with 4% paraformaldehyde (10 min)/permeabilized in 0.1% PBS-Tween used under the same conditions.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors