

Anti-Cytomegalovirus pp28 antibody [5C3] ab6502

12 References **2 图像**

概述

产品名称	Anti-Cytomegalovirus pp28抗体[5C3]
描述	小鼠单克隆抗体[5C3] to Cytomegalovirus pp28
宿主	Mouse
特异性	Cytomegalovirus pp28 tegument protein. Slight cross reactivity with HSV1.
经测试应用	适用于: WB, ELISA, ICC/IF
种属反应性	与反应: Human cytomegalovirus
免疫原	Tissue, cells or virus corresponding to Cytomegalovirus pp28.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: PBS
纯度	Protein A purified
纯化说明	> 95% IgG by SDS-PAGE.
克隆	单克隆
克隆编号	5C3
骨髓瘤	NS1/1-Ag4-1
同种型	IgG2a
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab6502于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/16000. Detects a band of approximately . kDa (predicted molecular weight: 28 kDa).
ELISA		1/102400.
ICC/IF		Use at an assay dependent concentration. PubMed: 20679731


靶标

相关性 Cytomegalovirus (CMV), is a genus of Herpes viruses; in humans the species is known as Human herpesvirus 5 (HHV-5). It belongs to the Betaherpesvirinae subfamily of Herpesviridae. The name means "very big cell virus".

图片



Please see protocols tab for procedure. Antigen is CMV infected cell extract coated at a dilution of 1/100.



See protocols tab for procedure. Antigen is CMV infected cell extract at 10 µg/cm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors