

Anti-DcR3 antibody [7G5] ab57956

[3 References](#) [2 图像](#)

概述

产品名称	Anti-DcR3抗体[7G5]
描述	小鼠单克隆抗体[7G5] to DcR3
宿主	Mouse
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human, Recombinant fragment
免疫原	Recombinant fragment, corresponding to amino acids 32-142 of Human DcR3
常规说明	<p>This product was changed from ascites to tissue culture supernatant on 12/3/19. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: PBS
纯度	Tissue culture supernatant
克隆	单克隆
克隆编号	7G5
同种型	IgG2a
轻链类型	kappa

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab57956于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration. This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.
IHC-P		Use at an assay dependent concentration.

靶标

功能	Decoy receptor for the cytotoxic ligands TNFS14/LIGHT and TNFSF6/FASL. Protects against apoptosis.
组织特异性	Detected in fetal lung, brain and liver. Detected in adult stomach, spinal cord, lymph node, trachea, spleen, colon and lung. Highly expressed in several primary tumors from colon, stomach, rectum, esophagus and in SW480 colon carcinoma cells.
序列相似性	Contains 4 TNFR-Cys repeats.
细胞定位	Secreted.


图片



DcR3 antibody (ab57956) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human tonsil.

This image was generated using the ascites version of the product.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DcR3 antibody [7G5] (ab57956)



Western blot against tagged recombinant protein immunogen using ab57956 DcR3 antibody at 1ug/ml. Predicted band size of immunogen is 38 kDa

This image was generated using the ascites version of the product.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors