

# Anti-DEFB-4 antibody ab70215

## 1 References

### 概述

产品名称	Anti-DEFB-4抗体
描述	兔多克隆抗体to DEFB-4
宿主	Rabbit
经测试应用	适用于: ICC/IF, Indirect ELISA, Sandwich ELISA, WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human DEFB-4. Highly pure (>98%). Database link: <a href="#">Q8WTQ1</a>
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Lyophilized:Centrifuge vial prior to opening. Reconstitute in sterile water to a concentration of 0.5 mg/ml. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab70215于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration. PubMed: 23627872
Indirect ELISA		Use at an assay dependent concentration. Can be paired for ELISA with Rabbit polyclonal to DEFB-4 (Biotin) ( <b>ab70216</b> ).
Sandwich ELISA		Use at an assay dependent concentration. Can be paired for Sandwich ELISA with <b>Biotin Rabbit polyclonal to DEFB-4 (ab70216)</b> . Capture
WB		Use at an assay dependent concentration.

## 靶标

功能	Has antimicrobial activity. Synergistic effects with lysozyme and DEFB103.
组织特异性	High expression in the testis. Gastric antrum exhibited relatively high levels. A lower expression is observed in uterus and neutrophils thyroid gland, lung, and kidney. No detectable expression in other tissues tested.
序列相似性	Belongs to the beta-defensin family.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors