

# Anti-Dnmt1 antibody [2C1] ab252825

重组

## 1 References [2 图像](#)


### 概述

产品名称	Anti-Dnmt1抗体[2C1]
描述	大鼠单克隆抗体[2C1] to Dnmt1
宿主	Rat
经测试应用	<b>适用于:</b> WB <b>不适用于:</b> IP
种属反应性	<b>与反应:</b> Human
免疫原	Recombinant fragment within Human Dnmt1 aa 1 to the C-terminus. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please <b>contact</b> our Scientific Support team to discuss your requirements. Database link: <a href="#">P26358</a>  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
阳性对照	WB: NCH1299, HEK-293, HeLa, Jurkat and HepG2 whole cell lysate.
常规说明	This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact <a href="mailto:orders@abcam.com">orders@abcam.com</a> . This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"><li>- High batch-to-batch consistency and reproducibility</li><li>- Improved sensitivity and specificity</li><li>- Long-term security of supply</li><li>- Animal-free production</li></ul> For more information <a href="#">see here</a> .

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: PBS, 0.05% BSA, 40% Glycerol (glycerin, glycerine)





**All lanes** : Anti-Dnmt1 antibody [2C1] (ab252825) at 0.353 µg/ml

**Lane 1** : NCI-H1299 (human lung carcinoma epithelial cell), whole cell lysate

**Lane 2** : HEK-293 (human embryonic kidney epithelial cell), whole cell lysate

**Lane 3** : HeLa (human cervix adenocarcinoma epithelial cell), whole cell lysate

**Lane 4** : Jurkat (human T cell leukemia T lymphocyte), whole cell lysate

**Lane 5** : HepG2 (human hepatocellular carcinoma epithelial cell), whole cell lysate

Lysates/proteins at 20 µg per lane.

**Secondary**

**All lanes** : Goat Anti-Rat IgG H&L (HRP) (**ab205720**) at 1/100000 dilution

**Predicted band size:** 183 kDa

**Exposure time:** 1 second

Blocking/Dilution buffer: 5% NFDM/TBST.

Why choose a recombinant antibody?

**Research with confidence**  
Consistent and reproducible results

**Long-term and scalable supply**  
Recombinant technology

**Success from the first experiment**  
Confirmed specificity

**Ethical standards compliant**  
Animal-free production

Anti-Dnmt1 antibody [2C1] (ab252825)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors