

Anti-DYNLT3 antibody ab121209

5 References **4 图像**

概述

产品名称	Anti-DYNLT3抗体
描述	兔多克隆抗体to DYNLT3
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Mouse, Rat, Human
免疫原	antigen sequence: DEVGFNAEEA HNVKECVDG VLGGEDYNHN NINQWTASV EQSLTHLVKL GKAYKYVTC AVVQKSAYGF HTASSCFWDT TSDGTCTVRW E, corresponding to amino acids 9-99 of Human DYNLT3. Run BLAST with Expasy Run BLAST with NCBI
阳性对照	WB: Wild-type HeLa, SK-MEL-20, MCF-7, and NBT-II whole cell lysates; IHC-P: Human cerebral cortex tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab121209于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/1000 - 1/2500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use a concentration of 0.04 - 0.4 µg/ml.

靶标

功能	Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. Probably binds BUB3 as part of transport cargo. Required for the efficient progression through mitosis.
序列相似性	Belongs to the dynein light chain Tctex-type family.
细胞定位	Nucleus. Cytoplasm > cytoskeleton. Chromosome > centromere > kinetochore. Colocalizes with BUB3 at kinetochores specifically during prometaphase.

图片



All lanes : Anti-DYNLT3 antibody (ab121209) at 0.04 µg/ml

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : DYNLT3 knockout HeLa cell lysate

Lane 3 : SH-SY5Y cell lysate

Lane 4 : U937 cell lysate

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

Observed band size: 13 kDa

False colour image of Western blot: Anti-DYNLT3 antibody (ab121209) staining at 0.04 µg/ml, shown in green; Mouse anti-Alpha Tubulin [DM1A] (ab7291) loading control staining at 1/20000 dilution, shown in red.

In Western blot, ab121209 was shown to bind specifically to

DYNLT3. A band was observed at 13 kDa in wild-type HeLa cell lysates with no signal observed at this size in DYNLT3 CRISPR-Cas9 edited cell lysates. The band observed in the CRISPR-Cas9 edited lysate lane above 13 kDa is likely to represent DYNLT3 with an insertion. This has not been investigated further and the functional properties of the gene product have not been determined. To generate this image, wild-type and DYNLT3 CRISPR-Cas9 edited HeLa cell lysates were analysed. First, samples were run on an SDS-PAGE gel then transferred onto a nitrocellulose membrane. Membranes were blocked in 3 % milk in TBS-0.1 % Tween® (TBS-T) before incubation with primary antibodies overnight at 4 °C. Blots were washed four times in TBS-T, incubated with secondary antibodies for 1 h at room temperature, washed again four times then imaged. Secondary antibodies used were Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed ([ab216776](#)) at 1/20000 dilution.




All lanes : Anti-DYNLT3 antibody (ab121209)

Lane 1 : SK-MEL-30 (human melanoma cell line) whole cell lysate


Lane 2 : MCF7 (human breast adenocarcinoma cell line) whole cell lysate

Loading control: Anti-PFN1.



ab121209 at 1/750 dilution staining DYNLT3 in Paraffin-embedded Human cerebral cortex tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DYNLT3 antibody (ab121209)



All lanes : Anti-DYNLT3 antibody (ab121209)

Lane 1 : NIH/3T3 (mouse embryo fibroblast cell line) whole cell lysate

Lane 2 : NBT-II (rat urinary bladder tumor cell line) whole cell lysate

Western blot - Anti-DYNLT3 antibody (ab121209)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors