

# Anti-FGF6 antibody [MM0290-8N21] ab89549

## 1 图像

### 概述

---

产品名称	Anti-FGF6抗体[MM0290-8N21]
描述	小鼠单克隆抗体[MM0290-8N21] to FGF6
宿主	Mouse
经测试应用	适用于: Neutralising, WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein (Human) FGF-6
阳性对照	WB: Human placenta tissue lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein G purified
纯化说明	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and filtered through a 0.2 µm filter.
克隆	单克隆
克隆编号	MM0290-8N21
同种型	IgG1

### 应用

---

## The Abpromise guarantee

**Abpromise™** 承诺保证使用ab89549于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Neutralising		Use at an assay dependent concentration.
WB		1/500 - 1/1000. Predicted molecular weight: 23 kDa.

## 靶标

功能	Plays an important role in the regulation of cell proliferation, cell differentiation, angiogenesis and myogenesis, and is required for normal muscle regeneration.
组织特异性	Leukemia cell lines with platelet/ megakaryocytic differentiation potential.
序列相似性	Belongs to the heparin-binding growth factors family.
细胞定位	Secreted > extracellular space.

## 图片



Anti-FGF6 antibody [MM0290-8N21] (ab89549) at 1/500 dilution + Human placenta tissue lysate

**Predicted band size:** 23 kDa

**Observed band size:** 27 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors