

FITC Mouse IgG2a, kappa monoclonal [MOPC-173] - isotype control ab18455

1 References

概述

产品名称	FITC荧光小鼠IgG2a, kappa单克隆抗体[MOPC-173]-同型对照
偶联物	FITC. Ex: 493nm, Em: 528nm
共轭说明	conjugated with FITC under optimal conditions.
经测试应用	适用于: Flow Cyt
免疫原	This information is considered to be commercially sensitive.
常规说明	<p>The ab18455, MOPC-173 immunoglobulin has unknown specificity and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed rat and human tissues.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Protein A purified
纯化说明	The ab18455, MOPC-173 antibody was purified by affinity chromatography The solution is free of unconjugated FITC.
Isotype control说明	The ab18455, MOPC-173 immunoglobulin has unknown specificity and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed rat and human tissues.
克隆	单克隆

克隆编号	MOPC-173
同种型	IgG2a
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab18455于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. Suitable for surface or intracellular staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors