

FITC-Rabbit IgG - Isotype Control ab37406

5 References 1 图像


概述	
产品名称	FITC荧光-Rabbit IgG -同型对照
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	This antibody has been selected to be an isotype control as it has no known specificity.
经测试应用	适用于: ELISA, Flow Cyt
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Store In the Dark.
存储溶液	pH: 7.40 Preservative: 0.1% Sodium azide Constituent: PBS
纯度	IgG fraction
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	Abpromise™ 承诺保证使用ab37406于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
Flow Cyt		Use 10µl for 10 ⁶ cells.

图片



FLISA plate was coated with Goat Anti-Rabbit IgG-UNLB and Goat Anti-Rabbit IgM-UNLB. Serially diluted Rabbit IgG-FITC ab37406 was captured and fluorescence intensity was quantified.

ELISA - FITC-Rabbit IgG - Isotype Control
(ab37406)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors