

Anti-Galectin 1 antibody [EPR3205] ab108389

重组 RabMAb

★★★★★ **1 Abreviews** **12 References** **6 图像**

概述

产品名称	Anti-Galectin 1抗体[EPR3205]
描述	兔单克隆抗体[EPR3205] to Galectin 1
宿主	Rabbit
经测试应用	适用于: WB, ICC/IF, IHC-P
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: HL60, K562, and NIH/3T3 cell lysates. Mouse heart, kidney, and spleen tissue lysates. Rat heart and kidney tissue lysates; ICC/IF: NIH/3T3 and K562 cells. IHC-P: Human placenta tissue.
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 9% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA, 50% Tissue culture supernatant</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR3205

同种型

IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab108389于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	1/1000 - 1/10000. Predicted molecular weight: 15 kDa.
ICC/IF		1/100.
IHC-P		1/3000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

靶标

功能

May regulate apoptosis, cell proliferation and cell differentiation. Binds beta-galactoside and a wide array of complex carbohydrates. Inhibits CD45 protein phosphatase activity and therefore the dephosphorylation of Lyn kinase.

组织特异性

Expressed in placenta, maternal decidua and fetal membranes. Within placenta, expressed in trophoblasts, stromal cells, villous endothelium, syncytiotrophoblast apical membrane and villous stroma. Within fetal membranes, expressed in amnion, chorioamniotic mesenchyma and chorion (at protein level). Expressed in cardiac, smooth, and skeletal muscle, neurons, thymus, kidney and hematopoietic cells.


序列相似性

Contains 1 galectin domain.

细胞定位


Secreted > extracellular space > extracellular matrix.

图片



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Galectin 1 antibody [EPR3205] (ab108389)

Immunohistochemical analysis of paraffin-embedded human placenta tissue labelling Galectin 1 with ab108389 at 1/3000 dilution (0.42 µg/ml) followed by a ready to use Rabbit specific IHC polymer detection kit HRP/DAB ([ab209101](#)) secondary antibody. Positive staining on vessels in human placenta tissue. The section was incubated with ab108389 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND® RX instrument. Counterstained with Hematoxylin. Heat mediated antigen retrieval with Citrate buffer (pH 6.0, epitope retrieval solution1) for 20 mins.



Western blot - Anti-Galectin 1 antibody [EPR3205] (ab108389)


All lanes : Anti-Galectin 1 antibody [EPR3205] (ab108389) at 1/1000 dilution

Lane 1 : HL60 cell lysate

Lane 2 : K562 cell lysate


Lysates/proteins at 10 µg per lane.

Predicted band size: 15 kDa




Immunocytochemistry/ Immunofluorescence - Anti-Galectin 1 antibody [EPR3205] (ab108389)

Immunocytochemistry/immunofluorescence analysis of K562 (human chronic myelogenous leukemia cell line from bone marrow) cells labelling Galectin 1 with purified ab108389 at 1/100. Cells were fixed with 4% paraformaldehyde and permeabilized with 0.1% Triton X-100. Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) ([ab150077](#)) at 1/1000 was used as the secondary antibody (green). Cells were counterstained with Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) ([ab195889](#)) at 1/200 dilution (red). Nuclear DNA was labelled with DAPI (blue).



Immunocytochemistry/ Immunofluorescence - Anti-Galectin 1 antibody [EPR3205] (ab108389)

Immunocytochemistry/immunofluorescence analysis of NIH/3T3 (mouse embryo fibroblast cell line) cells labelling Galectin 1 with purified ab108389 at 1/100. Cells were fixed with 4% paraformaldehyde and permeabilized with 0.1% Triton X-100. Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) ([ab150077](#)) at 1/1000 was used as the secondary antibody (green). Cells were counterstained with Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) ([ab195889](#)) at 1/200 dilution (red). Nuclear DNA was labelled with DAPI (blue).



Western blot - Anti-Galectin 1 antibody [EPR3205] (ab108389)

All lanes : Anti-Galectin 1 antibody [EPR3205] (ab108389) at 1/1000 dilution

Lane 1 : mouse heart tissue lysate

Lane 2 : mouse kidney tissue lysate

Lane 3 : mouse spleen tissue lysate

Lane 4 : rat heart tissue lysate

Lane 5 : rat kidney tissue lysate

Lane 6 : NIH 3T3 cell lysate

Lysates/proteins at 10 µg per lane.

Predicted band size: 15 kDa

Why choose a recombinant antibody?

Research with confidence
Consistent and reproducible results

Long-term and scalable supply
Recombinant technology

Success from the first experiment
Confirmed specificity

Ethical standards compliant
Animal-free production

Anti-Galectin 1 antibody [EPR3205] (ab108389)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors