

Anti-GFAP antibody [EPR1034Y] - Mouse IgG2a (Chimeric) - BSA and Azide free ab279302

重组

1 References 6 图像

概述

产品名称	Anti-GFAP抗体[EPR1034Y]-小鼠IgG2a (Chimeric) - BSA and Azide free
描述	小鼠单克隆抗体[EPR1034Y] to GFAP - Chimeric – BSA and Azide free
宿主	Mouse
经测试应用	适用于: WB, ICC, IHC-P, Flow Cyt (Intra), IHC-Fr, IP
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	IHC-P: Human cerebral cortex tissue. WB: Human, mouse and rat brain tissue lysate. IHC-Fr: Mouse cerebrum tissue. Flow Cyt (intra): Rat primary neural/glia cells. IP: Rat brain tissue lysate. ICC: mouse primary neural/glia cell
常规说明	<p>ab279302 is the carrier free version of ab279290.</p> <p>This mouse monoclonal chimeric antibody has been engineered from a RabMAb parent antibody (ab68428). By necessity, some rabbit sequence is retained as part of the variable domain. When multiplexing with other rabbit-derived antibodies, using cross absorbed Fc-reactive secondary antibodies are recommended.</p> <p>Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>This product is compatible with the Maxpar[®] Antibody Labeling Kit from Fluidigm, without the need for antibody preparation. Maxpar[®] is a trademark of Fluidigm Canada Inc.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Do Not Freeze.
存储溶液	pH: 7.2 Constituent: 100% PBS
无载体	是
纯度	Protein A purified
克隆	单克隆
克隆编号	EPR1034Y
同种型	IgG2a

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用 ab279302 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration.
ICC		1/50.
IHC-P		Use at an assay dependent concentration. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
Flow Cyt (Intra)		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration. Heat mediated antigen retrieval using sodium citrate buffer (10mM citrate pH 6.0 + 0.05% Tween-20).
IP		Use at an assay dependent concentration.

靶标

功能	GFAP, a class-III intermediate filament, is a cell-specific marker that, during the development of the central nervous system, distinguishes astrocytes from other glial cells.
组织特异性	Expressed in cells lacking fibronectin.
疾病相关	Defects in GFAP are a cause of Alexander disease (ALEXD) [MIM:203450]. Alexander disease is a rare disorder of the central nervous system. It is a progressive leukoencephalopathy whose hallmark is the widespread accumulation of Rosenthal fibers which are cytoplasmic inclusions in astrocytes. The most common form affects infants and young children, and is characterized by progressive failure of central myelination, usually leading to death usually within the first decade. Infants with Alexander disease develop a leukoencephalopathy with macrocephaly, seizures, and psychomotor retardation. Patients with juvenile or adult forms typically experience ataxia, bulbar signs and spasticity, and a more slowly progressive course.
序列相似性	Belongs to the intermediate filament family.


翻译后修饰

Phosphorylated by PKN1.

细胞定位

Cytoplasm. Associated with intermediate filaments.

图片




Immunocytochemistry - Anti-GFAP antibody
[EPR1034Y] - Mouse IgG2a (Chimeric) - BSA and
Azide free (ab279302)

This data was developed using **ab279290**, the same antibody clone in a different buffer formulation.

Immunocytochemical analysis of 4% PFA-fixed, 0.1% Triton X-100 permeabilized mouse neural/glia tissue labeling GFAP with **ab279290** at 1/50 (16 µg/ml) dilution followed by **ab98736** Goat Anti-Mouse IgG (DyLight® 488) pre-adsorbed at 1/1000 dilution (Green). Positive staining on mouse primary astrocytes is observed. The nuclear counterstain was DAPI (Blue).

Counterstain: **ab207165** anti-GFAP rabbit mAb at 1/1000 dilution with secondary antibody **ab150088** Goat Anti-Rabbit IgG H&L (Alexa Fluor® 594) pre-adsorbed at 1/1000 dilution.

Negative controls: **ab279290** with secondary **ab150080** at 1/500 dilution; and **ab207165** at 1/1000 dilution with secondary **ab98376** at 1/1000 dilution



Western blot - Anti-GFAP antibody [EPR1034Y] -
Mouse IgG2a (Chimeric) - BSA and Azide free
(ab279302)

All lanes : Anti-GFAP antibody [EPR1034Y] - Mouse IgG2a
(Chimeric) (**ab279290**) at 1/1000 dilution

Lane 1 : Human brain tissue lysate

Lane 2 : Mouse brain tissue lysate

Lane 3 : Rat brain tissue lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Peroxidase-Conjugated Goat anti-Mouse IgG (H+L) at
1/5000 dilution


Observed band size: 40-50 kDa

This data was produced using **ab279290**, the same clone in a different formulation.

Blocking/Dilution buffer: 5% NFDm/TBST.

The molecular weight observed is consistent with the literature (PMID: 824020, 2294, 6340792).

Exposure times: Lane 1: 3.25 seconds; Lane 2: 48 seconds; Lane 3: 10 seconds.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GFAP antibody [EPR1034Y] - Mouse IgG2a (Chimeric) - BSA and Azide free (ab279302)


This data was produced using [ab279290](#), the same clone in a different formulation.

IHC image of GFAP staining in a section of formalin-fixed paraffin-embedded normal human cerebral cortex performed on a Leica BOND™ system using the standard protocol F.

The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20mins. The section was then incubated with [ab279290](#), 1µg/ml, for 15 mins at room temperature. A rabbit anti-mouse IgG2a, was added for 8 mins at room temperature and detected using an HRP conjugated goat anti-rabbit compact polymer system. DAB was used as the chromogen.

The section was then counterstained with haematoxylin and mounted with DPX. The inset secondary-only control image is taken from an identical assay without primary antibody.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.




Immunohistochemistry (Frozen sections) - Anti-GFAP antibody [EPR1034Y] - Mouse IgG2a (Chimeric) - BSA and Azide free (ab279302)

This data was produced using [ab279290](#), the same clone in a different formulation.

Immunohistochemical analysis of 4% PFA-fixed, 0.2% Triton X-100 permeabilized frozen mouse cerebrum tissue labeling GFAP with [ab279290](#) at 1/500 (1.634 µg/ml) dilution followed by [ab150113](#) Goat Anti-Mouse IgG H&L (Alexa Fluor® 488) at 1/1000 dilution (Green). Positive staining on mouse cerebrum is observed. The nuclear counterstain was DAPI (Blue).

Secondary antibody control: Secondary antibody is [ab150113](#) Goat Anti-Mouse IgG H&L (Alexa Fluor® 488) at 1/1000 dilution.


Heat mediated antigen retrieval using sodium citrate buffer (10mM citrate pH 6.0 + 0.05% Tween-20).



Flow Cytometry (Intracellular) - Anti-GFAP antibody [EPR1034Y] - Mouse IgG2a (Chimeric) - BSA and Azide free (ab279302)

This data was produced using **ab279290**, the same clone in a different formulation.

Flow cytometric analysis of 4% paraformaldehyde fixed 90% methanol permeabilized rat primary neural/glia cells labelling GFAP with **ab279290** at 1/1000 dilution (0.1 µg)/ Right compared with a Mouse monoclonal IgG isotype control/ Left. Goat Anti-Mouse IgG (Alexa Fluor® 647, **ab150119**) at 1/2000 dilution was used as the secondary antibody.



Immunoprecipitation - Anti-GFAP antibody [EPR1034Y] - Mouse IgG2a (Chimeric) - BSA and Azide free (ab279302)

This data was produced using **ab279290**, the same clone in a different formulation.

GFAP was immunoprecipitated from 0.35 mg rat brain tissue lysate 10 µg with **ab279290** at 1/30 dilution (2 µg in 0.35mg lysates). Western blot was performed on the immunoprecipitate using **ab279290** at 1/1000 dilution. mouse IgG for IP (HRP) (**ab131368**) was used at 1/5000 dilution.

Lane 1: Rat brain tissue lysate 10 µg.

Lane 2: **ab279290** IP in rat brain tissue lysate.

Lane 3: Mouse monoclonal IgG2a (**ab18413**) instead of **ab279290** in rat brain tissue lysate.

Blocking and dilution buffer and concentration: 5% NFDm/TBST.

Exposure time: 23 seconds.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors