

Anti-GST antibody ab9085

★★★★★ [2 Abreviews](#) [59 References](#) [1 图像](#)

概述	
产品名称	Anti-GST抗体
描述	兔多克隆抗体to GST
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Species independent
免疫原	The details of the immunogen for this antibody are not available.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.00 Constituent: 50% Glycerol
纯度	Immunogen affinity purified
纯化说明	Affinity purified against bacterially expressed GST
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee **Abpromise™**承诺保证使用ab9085于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (2)	1/500.

靶标

相关性 GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

图片



All lanes : Anti-GST antibody (ab9085) at 1/2000 dilution

Lane 1 : Recombinant GST protein purified from E. coli XL1-Blue

Lane 2 : Recombinant GST-tagged p53 protein purified from E. coli XL1-Blue

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP-conjugated Goat anti-Rabbit IgG polyclonal at 1/20000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Observed band size: 26,80 kDa

Exposure time: 5 seconds

Blocked with 5% Milk for 1 hour at 25°C.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors